

# Steritest™ NEO Sterility Testing Devices

Convenience, Reliability and Safety

**Sterility testing is one of the most crucial steps in pharmaceutical product release.**

Our Steritest™ NEO devices simplify every aspect of testing, from handling to traceability, all within a closed concept system. The ease and convenience of this closed assembly enables you to increase productivity while maintaining the highest levels of quality and reliability. When used with the Steritest™ Symbio pump, the Steritest™ sterility test system delivers unmatched sterility testing consistency.

Since 1974, we are the market leader in sterility testing. The Steritest™ NEO devices, set a new standard for excellence, while maintaining all the advantageous membrane-base sealing technique.

## Benefits

- Thermo-sealed filtration membranes onto the base in all Steritest™ NEO units, ensure full integrity of the device, efficient membrane rinsing and eliminate the risk of by-pass
- Ergonomically designed needles fit the majority of test containers while maintaining a closed concept system
- Pre-installed colored clamps prevent any media filling errors and improve your workflow
- Canister design reduces foaming, enabling faster filtration
- Engraved information on each canister and peel-and-stick box label optimize traceability
- Volume graduation on the canisters improve your workflow accuracy (addition of a 25mL graduation)
- Pre-cut line on accessories bag to ease the opening
- Placement mark on tubing to ease the placement in the pump head\*

\* coming soon



### For Antibiotics and Products With Antimicrobial Agents

Ideally suited for antibiotic sample testing, the “Red Base” Steritest™ NEO devices incorporate our Durapore® (PVDF) membrane, offering broad chemical compatibility and low binding properties. The chemical composition of the extremely thin 0.45 µm Durapore® membrane provides low antibiotic binding, which improves rinsing efficiency and reduces the risk of false negatives. The canister design ensures efficient rinsing of residual antimicrobial agents and the needle and connections minimize the risk of antibiotic residuals.



#### Pre-installed Colored Clamps and Secure Tubing Length for Easy Handling

Steritest™ NEO canister tubing includes pre-installed colored clamps, eliminating the risk of error during media filling. Closing and re-opening the clamps during testing is as simple as pressing down or releasing as necessary. The 850 mm tubing length and the placement mark offer easier set-up with ample room and precision to perform the manipulations involved in sterility testing.

### For Products Without Antimicrobial Agents

The “Blue Base” Steritest™ NEO devices offer the ultimate in flexibility. Perfect for use with the majority of pharmaceutical drugs that do not have antimicrobial activity, our mixed esters of cellulose HA membrane allows fast flow rates for optimum throughput performance, and reduced testing time. A range of ergonomic needles is available to meet specific drug packaging requirements as well as simplify handling by gloved operators.




#### Canister Inlet Reduces Foaming and Speeds Filtration

In the past, analysts needed to reduce pump speed in order to minimize foaming when filtering protein solutions or adding media. The Steritest™ NEO canister inlet design reduces foaming, enabling faster filtration. This feature, coupled with our high-flow mixed esters of cellulose membrane, improves reliability while allowing you to maintain high pump speeds and reduce testing time.

### For Products requiring increased Chemical Compatibility

The “Green Base” Steritest™ NEO devices are suited for viscous products, such as creams and ointments, which are normally diluted in a sterile solvent, such as isopropyl myristate (IPM) to improve filterability. The solvent-resistant nylon canister and the Durapore® (PVDF) membrane ensure a perfect chemical compatibility with IPM and other solvents. The reinforced base structure and canister connection guarantee an unmatched resistance of the testing device during the whole filtration process. The “Green Base” Steritest™ NEO devices are the perfect choice for testing solvents, creams, ointments and veterinary injectables.



#### Needle Increase Efficiency and Security

Multiple needles adapters enable fast and easy manipulation of liquid drugs.

- Fits the majority of drug containers; reducing the risk of blockage when testing drugs in vials with thick septa
- Ergonomic, non-skid design provides a firm grip for gloved operations, especially in isolators

#### Closed Environment for Complete Testing Confidence

Using Steritest™ NEO devices ensures that pharmaceutical products are never exposed to the environment during the testing process. Sampling, filtering, rinsing, media transfer and incubating are all conducted within the Steritest™ NEO closed system. There is no need to open the containers or manipulate the membrane at any time—greatly reducing the risk of adventitious contamination.

#### Documented Qualification

We have compiled comprehensive Steritest Qualification Reports (available upon request) that confirms Steritest™ NEO device performance.

#### Consistent Performance

From the people who know sterility testing, we still take nothing for granted. We rigorously test each device during and after manufacturing.

- 100% integrity testing on every canister
- Strict physical and microbiological tests at every step of the assembly of the Steritest™ NEO device prior to release from manufacturing
- Certificate of quality provided with each system for your batch records
- Easy traceability with catalogue number, lot number, serial number and expiration date engraved on each canister

## Specifications

### Steritest™ NEO “Red Base” devices

For Antibiotics and Products Containing Antimicrobial Agents



|   |  |
|---|--|
| Canister Base Color   | Red  |
| Canister Base Membrane  | Low adsorption Durapore® membrane, 0.45 µm hydrophilic PVDF                            |
| Materials of Construction<br><i>Filtration Chamber (Canister):</i><br><i>Double Lumen Tubing:</i><br><i>Needle:</i> | Styrene acrylonitrile (SAN)<br>PVC, 850 mm length<br>Stainless steel and polyamide 6-6 |
| Sample Container Capacity   | 120 mL (graduation marks at 25, 50, 75 and 100 mL)                                     |
| Minimum Flow Rate (for water)   | 300 mL/min at 690 mbar (10 psi)  |
| Maximum Temperature   | 45 °C  |
| Maximum Operating Pressure  | 3.15 bars at 25 °C (45 psi at 77 °F)   |
| Sterilization   | Gamma irradiation  |

### Steritest™ NEO “Blue Base” devices

For Products without Antimicrobial Agents



|   |  |
|---|--|
| Canister Base Color   | Blue   |
| Canister Base Membrane  | Mixed Esters of Cellulose membrane, 0.45 µm                                      |
| Materials of Construction<br><i>Filtration Chamber (Canister):</i><br><i>Double Lumen Tubing:</i><br><i>Needle:</i> | Styrene acrylonitrile<br>PVC, 850 mm length<br>Stainless steel and polyamide 6-6 |
| Sample Container Capacity   | 120 mL (graduation marks at 25, 50, 75 and 100 mL)                               |
| Minimum Flow Rate (for water)   | 300 mL/min at 690 mbar (10 psi)  |
| Maximum Temperature   | 45 °C  |
| Maximum Operating Pressure  | 3.15 bars at 25 °C (45 psi at 77 °F)   |
| Sterilization   | Gamma irradiation  |

### Steritest™ NEO “Green Base” devices

For Products requiring increased chemical compatibility



|   |  |
|---|--|
| Canister Base Color   | Green  |
| Canister Base Membrane  | Low adsorption Durapore® membrane, 0.45 µm hydrophilic PVDF                      |
| Materials of Construction<br><i>Filtration Chamber (Canister):</i><br><i>Double Lumen Tubing:</i><br><i>Needle:</i> | polyamide 6-6 (nylon)<br>PVC, 850 mm length<br>Stainless steel and polyamide 6-6 |
| Sample Container Capacity   | 120 mL (graduation marks at 25, 50, 75 and 100 mL)                               |
| Minimum Flow Rate (for water)   | 300 mL/min at 690 mbar (10 psi)  |
| Maximum Temperature   | 45 °C  |
| Maximum Operating Pressure  | 3.15 bars at 25 °C (45 psi at 77 °F)   |
| Sterilization   | Gamma irradiation  |

## Ordering Information

### Steritest™ NEO devices for products without antimicrobial agents

| Description                                     | Pack size | Catalog number |
|---|-----------|----------------|
| <b>Liquids in ampoules and collapsible bags</b> |           |                |
| Steritest™ NEO device                           | 10        | TZHALA210      |
| Steritest™ NEO device double packed             | 10        | TZHALA205      |
| <b>Liquids in large vials</b>                   |           |                |
| Steritest™ NEO device                           | 10        | TZHALV210      |
| Steritest™ NEO device double packed             | 10        | TZHALV205      |
| <b>Liquids in small vials</b>                   |           |                |
| Steritest™ NEO device                           | 10        | TZHASV210      |
| Steritest™ NEO device double packed             | 10        | TZHASV205      |
| <b>Medical devices and collapsible bags</b>     |           |                |
| Steritest™ NEO device                           | 10        | TZHAMD210      |
| <b>Liquids in syringes</b>                      |           |                |
| Steritest™ NEO device                           | 10        | TZHASY210      |
| <b>Liquids in plastic containers</b>            |           |                |
| Steritest™ NEO device                           | 10        | TZHAPC210      |
| <b>Soluble powders in ampoules</b>              |           |                |
| Steritest™ NEO device                           | 10        | TZHADA210      |
| <b>Soluble powders in vials</b>                 |           |                |
| Steritest™ NEO device                           | 10        | TZHADV210      |

### Steridilutor® NEO devices and accessories

| Description   | Pack size | Catalog number |
|---|-----------|----------------|
| <b>Steridilutor® NEO device for vials</b>           |           |                |
| Steridilutor® NEO device without expansion chamber  | 10        | TZV000010      |
| Steridilutor® NEO device with expansion chamber     | 10        | TZVC00010      |
| <b>Steridilutor® NEO device for liquid transfer</b> |           |                |
| Steridilutor® NEO device for liquid transfer        | 10        | TZA000010      |
| <b>Sterile vent needles</b>                         |           |                |
| Sterile vent needles                                | 25        | TEFG02525      |

### Steritest™ NEO devices for antibiotics and products with antimicrobial agents

| Description                                     | Pack size | Catalog number |
|---|-----------|----------------|
| <b>Antibiotics</b>                              |           |                |
| Steritest™ NEO device                           | 10        | TZHVAB210      |
| Steritest™ NEO device double packed             | 10        | TZHVAB205      |
| <b>Powders and superpotent antibiotics</b>      |           |                |
| Steridilutor® NEO device with expansion chamber | 10        | TZVC00010      |
| Steritest™ NEO device                           | 10        | TZHVAB210      |
| <b>Recommended accessories</b>                  |           |                |
| Holder for Steridilutor® NEO vent chamber       | 1         | SYMBSVB01      |
| Sterile vent needles                            | 25        | TEFG02525      |
| <b>Liquids in large vial</b>                    |           |                |
| Steritest™ NEO device                           | 10        | TZHVLV210      |
| Steritest™ NEO device double packed             | 10        | TZHVLV205      |
| <b>Liquids in small vial</b>                    |           |                |
| Steritest™ NEO device                           | 10        | TZHVS210       |
| Steritest™ NEO device double packed             | 10        | TZHVS205       |
| <b>Soluble powders in vials</b>                 |           |                |
| Steritest™ NEO device                           | 10        | TZHVDV210      |
| Steritest™ NEO device double packed             | 10        | TZHVDV205      |
| <b>Medical devices and collapsible bags</b>     |           |                |
| Steritest™ NEO device                           | 10        | TZHVMD210      |

**Note:** For antibiotics or strong inhibitory products, the use of the TZHVAB210 filtration devices is highly recommended. Prior to the filtration step, it is recommended to dissolve and/or pool samples with the Steridilutor® NEO device.

### Steritest™ NEO devices for increased chemical compatibility

| Description  | Pack size | Catalog number |
|--|-----------|----------------|
| <b>Solvents, creams, ointments, and veterinary injectables</b> |           |                |
| Steritest™ NEO device  | 10        | TZHVSL210      |

Learn more at  
[EMDMillipore.com/steritestneo](https://EMDMillipore.com/steritestneo)

MilliporeSigma  
 400 Summit Drive  
 Burlington, MA 01803

Find contact information for your country at: [EMDMillipore.com/offices](https://EMDMillipore.com/offices)  
 For Technical Service, please visit: [EMDMillipore.com/techservice](https://EMDMillipore.com/techservice)

[EMDMillipore.com/biomonitoring](https://EMDMillipore.com/biomonitoring)

