

# Material Matters™

VOLUME 13 • NUMBER 2


## Atomic Layer Deposition Nanostructuring with ALD

NANOCOMPOSITE COATINGS with Tunable Properties Prepared by Atomic Layer Deposition

SILICON NITRIDE ATOMIC LAYER DEPOSITION: A Brief Review of Precursor Chemistry

GROUP 11 THIN FILMS by Atomic Layer Deposition

FEW MONOLAYER ATOMIC LAYER DEPOSITION on Surfaces and Interfaces for Energy Applications



## Introduction



**Taryn L. Fuhrman-Hall, M.B.A.**

Product Manager,  
Materials Science

Welcome to the second *Material Matters*™ issue of 2018. This issue explores atomic layer deposition (ALD) and the many disciplines in which it is utilized. The ALD process generates highly specific films with applications in many disciplines such as microelectronics, optics, sensors, solid-state detectors, battery and energy, nanomaterial growth, and solar power. The precision in which each layer can be deposited continues to make this a preferred process for advanced materials generation.

In the first article, Drs. Jeffrey Elam and Anil Mane (Argonne National Labs, USA) highlight the tunability of nanocomposite coatings and the remarkable suitability of atomic layer deposition for generating fit-for-purpose films. They specifically investigate the ALD mechanism of Mo:Al<sub>2</sub>O<sub>3</sub> and W:Al<sub>2</sub>O<sub>3</sub> films. They explore three specific practical applications of this technology for potential use in large-area photodetectors, electro-optical microsystems, and in concentrated solar power.

Silicon nitride is an essential material for semiconductor devices, and large-scale production of devices requires deposition of SiN at low temperatures. In the second article, Professor Jiyong Kim (University of Texas at Dallas, USA) reviews currently available silicon precursor classes: silanes, chlorosilanes, organosilanes, and heterosilanes. The relative film properties and growth characteristics for each class of precursor is discussed, with emphasis on new and up-and-coming chemistry.

The third article from Professor Sean Barry (Carleton University, Canada) reviews the accessibility of the three coinage metals: copper, silver, and gold. He reviews how the availability of precursors to deposit metal films is predicted by their standard potentials. The author discusses the well-known field of copper deposition, but also highlights the rapidly emerging area of silver and gold deposition.

Atomically precise modification of surfaces and interfaces with few monolayer material leads to significant enhancements in properties, performance and reliability of heterogeneous materials and devices. In the fourth article, Professor Parag Banerjee (Washington University, USA) highlights the impact of insulators based on a few monolayers and wide bandgap semiconductors and metals deposited using atomic layer deposition (ALD). These materials are used to boost the performance of photovoltaics, batteries and catalysts.

Each article in this publication is accompanied by a list of relevant materials available from MilliporeSigma Materials Science. For additional product information, visit us at [SigmaAldrich.com/matsci](https://www.sigmaaldrich.com/matsci). As always, please bother us with your new product suggestions as well as thoughts and comments for *Material Matters*™ at [matsi@sial.com](mailto:matsi@sial.com).

### About the Cover

Atomic layer deposition (ALD) uses alternating cycles of gaseous precursors to deposit materials onto a solid surface in an atomic layer-by-layer fashion. Our cover presents a literal and clear depiction of this layering effect in which the different color spheres represent the atoms of two distinct materials being layered onto a surface.

Merck KGaA  
Frankfurter Strasse 250  
64293 Darmstadt, Germany  
Phone +49 6151 72 0

#### To Place Orders / Customer Service

Contact your local office or visit  
[SigmaAldrich.com/order](https://www.sigmaaldrich.com/order)

#### Technical Service

Contact your local office or visit  
[SigmaAldrich.com/techinfo](https://www.sigmaaldrich.com/techinfo)

#### General Correspondence

Materials Science  
[materialsscience@sial.com](mailto:materialsscience@sial.com)

#### Subscriptions

Request your FREE subscription to *Material Matters*™ at [SigmaAldrich.com/mm](https://www.sigmaaldrich.com/mm)

The entire *Material Matters*™ archive is available at [SigmaAldrich.com/mm](https://www.sigmaaldrich.com/mm)

*Material Matters*™ (ISSN 1933-9631) is a publication of Merck KGaA.

Copyright (c) 2018 Merck KGaA, Darmstadt, Germany and/or its affiliates. All rights reserved. Merck, the vibrant M, Material Matters, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. More information on our branded products and services on [MerckMillipore.com](https://www.merckmillipore.com)

## Your Material Matters™



Bryce P. Nelson, Ph.D.  
Materials Science Initiative Lead


Dr. Alexander Spokoyny of the University of California — Los Angeles (USA) recommended the addition of cesium dodecahydroxydodecaborate (**902209**) to our catalog. This material represents a unique example of a robust inorganic polyol that can be used for the synthesis of well-defined small molecules, macromolecules and extended solid materials. It is conveniently synthesized and isolated as a Cs<sup>+</sup> salt from a reaction between B<sub>12</sub>H<sub>12</sub><sup>2-</sup> precursor and 30% H<sub>2</sub>O<sub>2</sub>.<sup>1-3</sup> Recently, B<sub>12</sub>(OH)<sub>12</sub><sup>2-</sup> has been used for the synthesis of weakly-coordinating photoredox active reagents capable of initiating olefin polymerization.<sup>4</sup> The stability of this cluster also allows the synthesis of atomically-precise nanomaterials capable of multivalent binding to protein surfaces<sup>5</sup> as well as extending solids featuring unique photochemical and electronic properties.<sup>6-7</sup>

## References

- (1) Peymann, T.; Knobler, C. B.; Khan, S. I.; Hawthorne, M. F. *J. Am. Chem. Soc.* **2001**, *123*, 2182–2185.
- (2) Bayer, M. J.; Hawthorne, M. F. *Inorg. Chem.* **2004**, *43*, 2018–2020.
- (3) Wixtrom, A. I.; Shao, Y.; Jung, D.; Machan, C. W.; Kevork, S. N.; Qian, E. A.; Khan, S. I.; Kubiak, C. P.; Spokoyny, A. M. *Inorg. Chem. Front.* **2016**, *3*, 711–717.
- (4) Messina, M. S.; Axtell, J. C.; Wang, Y.; Chong, P.; Wixtrom, A. I.; Kirlikovali, K. O.; Upton, B. M.; Hunter, B. M.; Shafaat, O. S.; Khan, S. I.; Winkler, J. R.; Gray, H. B.; Alexandrova, A. N.; Maynard, H. D.; Spokoyny, A. M. *J. Am. Chem. Soc.* **2016**, *138*, 6952–6955.
- (5) Qian, E. Q.; Wixtrom, A. I.; Axtell, J. C.; Saebi, A.; Rehak, P.; Han, Y.; Mouilly, E. H.; Mosallaei, D.; Chow, S.; Messina, M.; Wang, J.-Y.; Royappa, A. T.; Rheingold, A. L.; Maynard, H. D.; Kral, P.; Spokoyny, A. M. *Nature Chem.* **2017**, *9*, 333–340.
- (6) Jung, D.; Saleh, L. M. A.; Berkson, Z.; El-Kady, M. F.; Hwang, J. Y.; Mohamed, N.; Wixtrom, A. I.; Titarenko, E.; Shao, Y.; McCarthy, K.; Guo, J.; Martini, I. B.; Kraemer, S.; Wegener, E. C.; Saint-Cricq, P.; Ruehle, B.; Langeslay, R. R.; Delferro, M.; Brosmer, J. L.; Hendon, C. H.; Gallagher-Jones, M.; Rodriguez, J.; Chapman, K. W.; Miller, J. T.; Duan, X.; Kaner, R. B.; Zink, J. I.; Chmelka, B. F.; Spokoyny, A. M. *Nat. Mater.* **2018**, *17*, 341–348.
- (7) Axtell, J. C.; Saleh, L. M. A.; Qian, E. A.; Wixtrom, A. I.; Spokoyny, A. M. *Inorg. Chem.* **2018**, *57*, 2333–2350.

## Cesium dodecahydroxydodecaborate

[12589-26-3] B<sub>12</sub>H<sub>12</sub>O<sub>12</sub>·2Cs



Boric acid (H<sub>14</sub>B<sub>12</sub>O<sub>12</sub>), dicesium salt; Dodecaborate(2-), dodecahydroxy-, dicesium (9Cl); Dicesium dodecahydroxy-closo-dodecaborate; Dicesium dodecahydroxydodecaborate(2-).

Store at 2-8°C

1g

902209-1G

## Table of Contents

## Articles

Nanocomposite Coatings with Tunable Properties Prepared by Atomic Layer Deposition	45
Silicon Nitride Atomic Layer Deposition: A Brief Review of Precursor Chemistry	55
Group 11 Thin Films by Atomic Layer Deposition	60
Few Monolayer Atomic Layer Deposition on Surfaces and Interfaces for Energy Applications	65

## Featured Products

Precursors Packaged for Deposition Systems A selection of available ALD materials	53
Sputtering Targets A list of targets	54
Silicon Precursors A selection of ALD and solution deposition materials	57
Group 11 Deposition Materials A list of copper, silver and gold deposition materials	64
Titania Materials for Support and Solar Cells A list of titania nanomaterials and precursors	70
Lithium Precursors for Energy A selection of lithium precursors	70
Zinc Precursors for Catalysis A selection of zinc precursors	71

# Beyond pure



Ultra-high purity metals, salts, and oxides are essential in creating advanced materials. Explore the breadth of our portfolio's materials with purity of 99.999% or higher for applications such as:

- Photovoltaics
- Phosphor materials
- Nanoparticle synthesis
- Magnetic memory
- Biomedical applications
- Electronic devices

Product Description	Cat. No.
<b>Metals</b>	
Aluminum, wire, diam. 1.0 mm, 99.999% trace metals basis	266558
Antimony, beads, 1-5 mm, low oxide, 99.999% trace metals basis	452343
Bismuth, pieces, 1-12 mm, 99.999% trace metals basis	556130
Copper, powder, 99.999% trace metals basis	203122
Gallium, 99.9995% trace metals basis	203319
Germanium, chips, 99.999% trace metals basis	203343
Gold, beads, 1-6 mm, 99.999% trace metals basis	326542
Indium, beads, diam. 2-5 mm, 99.999% trace metals basis	264113
Selenium, pellets, <5 mm particle size, ≥99.999% trace metals basis	204307
Silicon, powder, -60 mesh, 99.999% trace metals basis	267414
Zinc, foil, thickness 0.25 mm, 99.999% trace metals basis	267619
<b>Oxides</b>	
Bismuth(III) oxide, powder, 99.999% trace metals basis	202827
Copper(II) oxide, 99.999% trace metals basis	203130
Europium(III) oxide, 99.999% trace metals basis	323543
Germanium(IV) oxide, powder, 99.999% trace metals basis	483001
Holmium(III) oxide, powder, 99.999% trace metals basis	229679
Iron(III) oxide, ≥99.995% trace metals basis	529311
Lanthanum(III) oxide, 99.999% trace metals basis	203556
Lead(II) oxide, 99.999% trace metals basis	203610
Mercury(II) oxide, 99.999% trace metals basis	203793
Zinc oxide, 99.999% trace metals basis	204951

Product Description	Cat. No.
<b>Salts</b>	
Aluminum chloride, anhydrous, powder, 99.999% trace metals basis	563919
Ammonium acetate, 99.999% trace metals basis	372331
Barium chloride dihydrate, ≥99.999% trace metals basis	529591
Cesium chloride, ≥99.999% trace metals basis	203025
Copper(I) bromide, 99.999% trace metals basis	254185
Copper(II) nitrate hydrate, 99.999% trace metals basis	229636
Gallium(III) chloride, beads, anhydrous, ≥99.999% trace metals basis	427128
Gold(III) chloride hydrate, 99.995% trace metals basis	254169
Indium(III) chloride, 99.999% trace metals basis	203440
Lanthanum(III) nitrate hexahydrate, 99.999% trace metals basis	203548
Lead(II) bromide, 99.999% trace metals basis	398853
Lead(II) iodide, beads, -10 mesh, 99.999% trace metals basis	554359
Silver nitrate, 99.9999% trace metals basis	204390
Sodium chloride, 99.999% trace metals basis	204439
Sodium nitrite, 99.999% trace metals basis	563218
Zinc nitrate hydrate, 99.999% trace metals basis	230006

For a complete list of ultra-high purity products, visit [SigmaAldrich.com/uhp](http://SigmaAldrich.com/uhp)



# Nanocomposite Coatings with Tunable Properties Prepared by Atomic Layer Deposition



Anil U. Mane and Jeffrey W. Elam\*

Argonne National Laboratory, Argonne, Illinois 60439, USA

\*Email: jelam@anl.gov

## Introduction

Nanocomposite coatings, comprised of a precisely blended mixture of components, have diverse applications in microelectronics, optics, sensors, and solid-state detectors. By varying the composition of the mixture, the electrical, optical, and physical properties of the coatings can be “tuned” over almost the full range of the individual components, and can sometimes yield properties distinct from those of any of the constituents. There are many techniques for producing nanocomposite coatings, but atomic layer deposition (ALD) is particularly well suited for controlling the properties of the fabricated coatings. ALD uses alternating cycles of gaseous precursors to a solid surface to deposit materials in an atomic layer-by-layer fashion.<sup>1</sup> Because a chemical reaction on the surface naturally terminates when the surface functional groups are all reacted, the amount of material deposited in one ALD cycle is easily controlled. This self-limiting property, coupled with facile diffusion of the precursor vapors into narrow pores and voids, allows complex, 3-dimensional substrates to be coated with excellent uniformity and conformality.<sup>2</sup> ALD can be used to deposit a wide variety of materials including metal oxides, nitrides, sulfides, and even pure elements. By alternating between two materials (e.g., a metal and a metal oxide), we can synthesize complex nanocomposite coatings.<sup>3-5</sup> The thickness of the coating is controlled by the total number of ALD cycles performed, and the composition is controlled by the ratio of the ALD cycles executed for each of the two components.

In this study, we employed *in situ* quartz crystal microbalance (QCM) measurements and Fourier transform infrared (FTIR) absorption spectroscopy to investigate the ALD growth mechanism of nanocomposite films comprised of Mo:Al<sub>2</sub>O<sub>3</sub> and W:Al<sub>2</sub>O<sub>3</sub>. These films were prepared using alternating exposures to trimethyl aluminum (TMA) and H<sub>2</sub>O for the Al<sub>2</sub>O<sub>3</sub> ALD, and

alternating MF<sub>6</sub>/Si<sub>2</sub>H<sub>6</sub> exposures for the metal ALD, where M=Mo or W. In the case of the Mo:Al<sub>2</sub>O<sub>3</sub> films, QCM showed that the Mo ALD inhibits the Al<sub>2</sub>O<sub>3</sub> ALD and vice versa. Despite this inhibition, the relationship between Mo content and Mo cycle percentage was easily controlled. Surprisingly, FTIR revealed that the reducing agent for the Mo was not the Si<sub>2</sub>H<sub>6</sub>, but rather the TMA exposure from the subsequent Al<sub>2</sub>O<sub>3</sub> ALD cycle. Elemental analysis showed that the M:Al<sub>2</sub>O<sub>3</sub> films were uniform in composition and contained Al, O, and metallic Mo or W as expected, but also included significant F and C. Cross-sectional transmission electron microscopy revealed the film structure to be metallic nanoparticles (~1 nm) embedded in an amorphous matrix.

Three distinct practical applications of this technology are explored in this paper. We utilized these nanocomposite coatings to functionalize capillary glass array plates to fabricate large-area microchannel plates (MCPs) suitable for application in large-area photodetectors. In addition, we applied these films as charge drain coatings in electron-optical microsystems for a prototype electron beam lithography tool, permitting high-resolution electron beam patterns without charging artifacts. Finally, we explored the suitability of these films for use as selective solar absorbing coatings in concentrated solar power.

## Experimental

Deposition of the Mo:Al<sub>2</sub>O<sub>3</sub> and W:Al<sub>2</sub>O<sub>3</sub> composite ALD film was performed at 200 °C in a hot-walled viscous flow ALD reactor.<sup>6</sup> This reactor was equipped with a quartz crystal microbalance (QCM) that allowed *in situ* studies of the ALD processes. The TMA (97%, Sigma-Aldrich), deionized H<sub>2</sub>O, Si<sub>2</sub>H<sub>6</sub> (99.998%, Sigma-Aldrich), MoF<sub>6</sub> (99%, Alfa Aesar), and WF<sub>6</sub> (99.8%, Sigma-Aldrich) precursors were maintained at room temperature. The

ultrahigh purity (99.999%)  $N_2$  carrier gas flow was set to 300 sccm, which provided a base pressure of 1.0 Torr in the ALD reaction chamber as measured by a heated Baratron pressure gauge (MKS model 629B). The  $Mo:Al_2O_3$  and  $W:Al_2O_3$  composite films were deposited on  $n$ -type Si(100) substrates. Prior to ALD, the substrates were cleaned using a 10 min ultrasonic treatment in acetone. For the  $Al_2O_3$  ALD, TMA and  $H_2O$  were alternately pulsed into the  $N_2$  carrier flow with the following sequence: 1s TMA dose (0.2 Torr) – 5s purge – 1s  $H_2O$  dose (0.3 Torr) – 5s purge. Similarly, the Mo and W ALD used alternating exposures to  $MoF_6/Si_2H_6$  and  $WF_6/Si_2H_6$ , respectively, with the timing sequence: 1s  $Si_2H_6$  dose (0.25 Torr) – 5s purge – 1s  $MoF_6$  or  $WF_6$  dose (0.05 Torr) – 5s purge. These conditions for the  $Al_2O_3$ , Mo, and W ALD provided self-limiting growth as verified by *in situ* QCM measurements.

*In situ* QCM was performed to examine the  $Mo:Al_2O_3$  and  $W:Al_2O_3$  composite ALD using different ALD pulse sequences. The QCM measurements typically used 10s  $N_2$  purge times to allow the QCM signal to stabilize after each precursor exposure. Transmission electron microscopy (TEM) analysis was performed by Evans Analytical Group (Sunnyvale, CA). TEM samples were prepared using the *in situ* focussed ion beam (FIB) lift out technique on an FEI Strata Dual Beam FIB/SEM. The samples were capped with a protective layer of carbon prior to FIB milling, and were imaged with a FEI Tecnai TF-20 FEG/TEM operated at 200 kV in bright-field (BF) TEM mode, high-resolution (HR) TEM mode, and nanobeam diffraction (NBD) mode. The composition of the  $Mo:Al_2O_3$  and  $W:Al_2O_3$  composite layers was determined by depth profiling using X-ray photoelectron spectroscopy (XPS, Evans Analytical Group) and Rutherford backscattering spectroscopy (RBS, Evans Analytical Group). The metal content in the composite films was measured using X-ray fluorescence (XRF, Oxford ED2000).

The resistance of the  $Mo:Al_2O_3$  and  $W:Al_2O_3$  composite layers was determined by performing current-voltage (I-V) measurements using a Keithley Model 6487 current-voltage source. To facilitate I-V measurements on these highly resistive coatings, the films were deposited on insulating substrates with lithographically patterned comb structures comprised of interdigitated Au electrodes with a 2  $\mu m$  spacing.<sup>7</sup> These comb structures provided an 80,000 times greater effective contact area, and an equivalent boost in current compared to conventional four-point probe measurements. Additional I-V and electron amplification measurements were performed on  $W:Al_2O_3$  layers deposited on high aspect ratio (60:1) borosilicate glass capillary arrays with 20  $\mu m$  pores fabricated by Incom, Inc. (Charlton, MA). These capillary arrays were also coated with an eight nm ALD MgO emissive layer to boost the secondary electron coefficient, and both sides were coated with 100 nm evaporated nickel-chromium to provide a low resistance electrical contact.


*In situ* Fourier transform infrared (FTIR) measurements were performed during the  $Mo:Al_2O_3$  ALD using a Nicolet 6700 with a liquid nitrogen-cooled MCT-B detector. For these measurements, the FTIR beam was diverted through CsI windows and into the

ALD reactor. The windows were protected by gate valves, which were closed during the ALD exposures to prevent deposition on the windows. The FTIR beam was directed through a sample comprised of  $ZrO_2$  nanopowder pressed into a stainless steel grid to enhance the sample surface area and boost the optical absorption.<sup>8</sup> Prior to the  $Mo:Al_2O_3$  ALD, the  $ZrO_2$  powder was coated using 10 cycles of TMA and  $H_2O$  for  $Al_2O_3$  ALD.

## Results and Discussion


### Growth and Properties of $Mo:Al_2O_3$ Composite Films

Figure 1 shows a schematic of the precursor dosing strategy used to synthesize the  $Mo:Al_2O_3$  ALD composite films. Alternating TMA/ $H_2O$  exposures were used to deposit  $Al_2O_3$ , and periodically a single  $MoF_6/Si_2H_6$  cycle was performed to introduce the Mo. We define the Mo cycle percentage as:  $\%Mo = Mo / (Mo + Al_2O_3) * 100$  where Mo and  $Al_2O_3$  are the relative numbers of TMA/ $H_2O$  and  $MoF_6/Si_2H_6$  cycles performed, respectively. Because of the drastically different resistivity for these two materials ( $\sim 10^{16}$  and  $10^{-4}$  Ohm cm for  $Al_2O_3$  and Mo, respectively), we expected to be able to tune the resistivity for the composite film over a broad range by adjusting the Mo cycle percentage.




**Figure 1.** Schematic illustration of precursor pulsing sequence for ALD of  $Mo:Al_2O_3$  nanocomposite coating using alternating TMA/ $H_2O$  exposures for the  $Al_2O_3$  ALD and alternating  $MoF_6/Si_2H_6$  exposures for the Mo ALD.

Figure 2A shows *in situ* QCM measurements performed during the ALD of a  $Mo:Al_2O_3$  composite film using a 10% Mo cycle percentage (nine  $Al_2O_3$  ALD cycles followed by one Mo ALD cycle). The overall pattern consists of nine small mass gains followed by a single, larger mass gain, and this pattern is repeated three times. Figure 2B shows an expanded view of the boxed region in Figure 2A and reveals some informative details about the QCM mass gains. Between 95–115s, two TMA/ $H_2O$  cycles for  $Al_2O_3$  ALD are shown and the QCM signals match the expected pattern for  $Al_2O_3$  ALD.<sup>3</sup> The net mass gain for each  $Al_2O_3$  ALD cycle is  $\sim 30$  ng/cm<sup>2</sup>, which is slightly smaller than the expected value. In contrast, the QCM signals observed during the  $MoF_6/Si_2H_6$  exposures do not follow the pattern reported in the literature for Mo ALD.<sup>9</sup> In particular, no mass gain is observed during the  $Si_2H_6$  exposure, and the net mass gain is only  $\sim 250$  ng/cm<sup>2</sup>, which is far below the 1,000 ng/cm<sup>2</sup> expected for Mo ALD. It is interesting to note that the first TMA exposure following the Mo ALD cycle yields a sharp mass rise followed by a slow decay. This is the signature pattern for a temperature-induced QCM transient produced by an exothermic reaction.<sup>10</sup>




**Figure 2.** *In situ* QCM measurements during Mo:Al<sub>2</sub>O<sub>3</sub> composite film ALD using a 10% Mo cycle percentage (nine Al<sub>2</sub>O<sub>3</sub> ALD cycles followed by one Mo ALD cycle). **A)** Mass vs. time over 30 ALD cycles. **B)** Expanded view of the inset in **A**, with TMA, H<sub>2</sub>O, Si<sub>2</sub>H<sub>6</sub>, and MoF<sub>6</sub> dose times are indicated.

**Figure 3** shows the net mass gains following each ALD cycle for Mo (triangles) and Al<sub>2</sub>O<sub>3</sub> (circles). The dashed, horizontal line shows the expected, steady-state mass gains for Al<sub>2</sub>O<sub>3</sub> ALD at ~35 ng/cm<sup>2</sup>. This figure illustrates that the Al<sub>2</sub>O<sub>3</sub> ALD is inhibited by ~50% initially following the Mo ALD, and does not quite achieve the steady-state value even after nine Al<sub>2</sub>O<sub>3</sub> ALD cycles. Moreover, the Mo ALD is also severely inhibited on the Al<sub>2</sub>O<sub>3</sub> surface by ~75%.




**Figure 3.** Net mass gains following each ALD cycle of Mo:Al<sub>2</sub>O<sub>3</sub> composite film deposition using a 10% Mo cycle percentage. Mo and Al<sub>2</sub>O<sub>3</sub> ALD cycles are indicated by the triangles and circles, respectively. The dashed, horizontal line shows the expected steady-state mass gain of 35 ng/cm<sup>2</sup> for Al<sub>2</sub>O<sub>3</sub> ALD.

The origin for some of the unusual QCM results can be gleaned from the *in situ* FTIR measurements. **Figure 4** shows the IR absorbance at 4000 cm<sup>-1</sup>, measured following the individual precursor exposures, where the precursor pulsed immediately prior to recording the IR spectrum is indicated above each data point. The 4000 cm<sup>-1</sup> frequency was selected because it falls in a spectral region where there are no peaks from any of the surface functional groups, and therefore serves as an indirect measure of the relative conductivity of the film.<sup>11</sup> The IR absorbance at 4000 cm<sup>-1</sup> remains relatively constant during the TMA and H<sub>2</sub>O exposures for Al<sub>2</sub>O<sub>3</sub> ALD. We might expect the IR absorbance to increase during the Si<sub>2</sub>H<sub>6</sub> exposures for the Mo ALD since Si<sub>2</sub>H<sub>6</sub> is believed to be the reducing agent for the Mo.<sup>9</sup> However, this is clearly not the case, as there is practically no change in the IR absorbance during the Si<sub>2</sub>H<sub>6</sub> exposures. Instead, the IR absorbance increases sharply during the TMA exposures that follow the MoF<sub>6</sub> exposures. This finding suggests that TMA is in fact the reducing agent, which is consistent with the exothermic transients seen in the QCM data during the TMA exposures following MoF<sub>6</sub> exposures (**Figure 2B**). Additional evidence for the reducing effect of TMA on the adsorbed MoF<sub>6</sub> was indicated by a decrease in the MoOF<sub>x</sub> stretching frequency at 1040 cm<sup>-1</sup> during the TMA exposures (spectra not shown).



**Figure 4.** IR absorbance at 4000 cm<sup>-1</sup> measured by *in situ* FTIR following individual precursor exposures for the Mo:Al<sub>2</sub>O<sub>3</sub> composite film ALD. The precursor pulsed immediately prior to recording the IR spectrum is indicated above each data point.


Figure 5 shows a low-resolution cross-sectional TEM image of a 40 nm Mo:Al<sub>2</sub>O<sub>3</sub> composite film prepared on a Si(100) substrate using 10% Mo cycles. This film appears dense and continuous, and the top surface is relatively smooth and parallel to the Si substrate. The interface between the Si substrate and the Mo:Al<sub>2</sub>O<sub>3</sub> composite film shows a low density, amorphous region attributed to the Si native oxide and Al<sub>2</sub>O<sub>3</sub> from the initial ALD cycles. Under higher resolution (Figure 5A), the Mo:Al<sub>2</sub>O<sub>3</sub> composite film is seen to consist of 1-2 nm particles (dark spots) embedded in a lower density matrix. Close inspection of the particles (Figure 5A, inset) reveals weak lattice fringes. Moreover, nano-beam diffraction measurements acquired from the nanocomposite film (Figure 5B) exhibit diffuse rings consistent with crystalline nanoparticles. Based on these observations, and on our previous study of the W:Al<sub>2</sub>O<sub>3</sub> composite film ALD,<sup>5</sup> we hypothesize that the crystalline nanoparticles in Figure 5A are metallic Mo. These Mo nanoparticles might form through the sintering of Mo atoms upon reduction of MoOF<sub>x</sub> surface species by the TMA. The role of TMA as a reducing agent is somewhat surprising, but we believe this to be true based on our FTIR results, and because the microstructure and properties of the Mo:Al<sub>2</sub>O<sub>3</sub> composite films are practically identical when the Si<sub>2</sub>H<sub>6</sub> exposures are not performed.



**Figure 5.** Cross-sectional TEM analysis of 40 nm Mo:Al<sub>2</sub>O<sub>3</sub> composite film prepared on Si(100) substrate using 10% Mo cycles. **A)** High resolution image showing 1-2 nm particles (dark spots) embedded in a lower density matrix. **B)** Nano-beam diffraction of film region shown in **A**. **C)** Low resolution image of entire film thickness.

The composition of the Mo:Al<sub>2</sub>O<sub>3</sub> composite films was determined using XPS depth profiling and RBS measurements of films prepared on Si(100) substrates. These measurements revealed that the films contained Al, O, and Mo as expected, but also contained ~5% C and ~10% F. XPS demonstrated that Mo exists in both the metallic and sub-oxide (Mo<sup>+4</sup>) states in the film. These findings, when combined with the QCM and FTIR measurements, suggest that the TMA reduces the adsorbed MoOF<sub>x</sub> species to form AlF<sub>3</sub> and metallic Mo. In addition, the C and F presence in the films is similar to that in a previous study which concluded that alternating exposures to NbF<sub>5</sub> and TMA produced films comprised of NbF<sub>x</sub>, NbC, C, and AlF<sub>3</sub>.<sup>12</sup>

Next, a series of Mo:Al<sub>2</sub>O<sub>3</sub> composite films was prepared on Si(100) substrates using different Mo cycle percentages and analyzed by X-ray fluorescence to determine the Mo content of the films. As shown in Figure 6A, the Mo content increases linearly from 2.1 to 9.0 mol% Mo as the Mo cycle percentage increases from eight to 15% Mo. This linear relationship demonstrates that the Mo content can be easily and precisely controlled by adjusting the Mo cycle percentage during the Mo:Al<sub>2</sub>O<sub>3</sub> composite film ALD. The resistivity of the Mo:Al<sub>2</sub>O<sub>3</sub> composite films was determined using current-voltage measurements of films deposited on comb structures, and the results are shown in Figure 6B. Figure 6B shows that the resistivity of the Mo:Al<sub>2</sub>O<sub>3</sub> composite films can be tuned over a very broad range by adjusting the Mo cycle percentage. The range of resistivity values provided by these coatings (10<sup>4</sup>–10<sup>10</sup> Ohm cm) is particularly useful for applications in electron multipliers and charge drain coatings, and this capability will be demonstrated below.




**Figure 6.** **A)** Mo content of ALD Mo:Al<sub>2</sub>O<sub>3</sub> composite films as determined by XRF measurements versus Mo cycle percentage. **B)** Resistivity of ALD Mo:Al<sub>2</sub>O<sub>3</sub> composite films as determined by current-voltage measurements versus Mo cycle percentage.



## Growth and Properties of W:Al<sub>2</sub>O<sub>3</sub> Composite Films


We performed *in situ* QCM measurements of the W:Al<sub>2</sub>O<sub>3</sub> composite film grown by ALD to explore their growth behavior. We deposited and characterized a range of W:Al<sub>2</sub>O<sub>3</sub> composite films while varying the W cycle percentage to determine the thickness, composition, morphology, and microstructure of these films and to establish the effect of W content on resistivity.<sup>5</sup> Some of these results are presented in Figure 7.



**Figure 7.** A) W content of ALD W:Al<sub>2</sub>O<sub>3</sub> composite films deduced from *in situ* QCM versus W cycle percentage. B) Resistivity of ALD W:Al<sub>2</sub>O<sub>3</sub> composite films determined by current-voltage measurements versus W cycle percentage.

Figure 7A shows the composition of the W:Al<sub>2</sub>O<sub>3</sub> composite films as determined by *in situ* QCM versus W cycle percentage. The W mol% varies smoothly over the full range, and is approximately linear between 10–30%. Figure 7B shows the effect of varying the W cycle percentage on the resistivity as determined by current-voltage measurements of films deposited on comb chips and on capillary glass array substrates. The resistivity decreases approximately exponentially from  $\sim 10^{12}$  Ohm cm at 10%W cycles to  $\sim 10^8$  Ohm cm at 30% W cycles. This behavior is similar to that of the Mo:Al<sub>2</sub>O<sub>3</sub> composite films (Figure 6B), although the slope is much steeper in the case of the Mo:Al<sub>2</sub>O<sub>3</sub> composite films (6 decades over 4% change in Mo cycles) as compared to W:Al<sub>2</sub>O<sub>3</sub> composite films (4 decades over 20% change in W cycles). This difference may stem from the much higher deposition rate for the Mo ( $\sim 10$  Å/cycle)<sup>9</sup> versus the W ( $\sim 4$  Å/cycle)<sup>13</sup> under similar growth conditions, especially in light of the exponential relationship between resistivity and metal content.

Next, we performed optical absorption measurements on the W:Al<sub>2</sub>O<sub>3</sub> composite films prepared using 0–60% W cycles. The direct bandgaps were extracted from the absorption spectra using Tauc plots, and the results are shown in Figure 8A. The bandgap decreases with increasing W cycle percentage. Figure 8B shows the absorption coefficients as a function of photon energy for these coatings and clearly shows that the absorption onset shifts to lower energies with increasing W cycle percentages between 10–40%, and then transitions to a broadband, metallic absorption for the higher W cycle percentages of 50–60%. This tunable absorption property is desirable for selective solar absorption coatings that must absorb visible light while transmitting or reflecting infrared wavelengths, as discussed next.



**Figure 8.** A) Direct bandgap versus W cycle percentage and B) absorption coefficient versus photon energy for ALD W:Al<sub>2</sub>O<sub>3</sub> composite films prepared using 0–60% W cycle percentage.

## Other Nanocomposite Coatings

Similar to the coating strategies previously described, we have developed ALD methods for nanocomposite coatings using combinations of W, Mo, Ta, Nb, and Co metals, and Al<sub>2</sub>O<sub>3</sub>, ZrO<sub>2</sub>, HfO<sub>2</sub>, MgO, and TiO<sub>2</sub> dielectrics. Although a complete description is beyond the scope of this article, we note that the trends of increasing electrical conductivity, decreasing bandgap, and increasing absorption coefficient with increasing metal content are common to all of these materials.

## Applications for Tunable Resistance Coatings


### Microchannel Plates (MCPs)

Microchannel plates (MCPs) are solid-state electron amplifiers comprised of a two-dimensional array of microscopic channels in the form of a thin, flat plate. MCPs are used in detectors to identify low levels of electrons, ions, photons, or neutrons, and provide an amplified response via gain from secondary electron emission that occurs within the individual channels of the MCP.<sup>14</sup> This unique quality allows MCPs to form images, and is critical to a wide variety of detectors and devices. Conventional MCP fabrication involves using multi-fiber glass working techniques to draw and assemble an array of solid-core fibers into a block. The block is subsequently sliced into thin wafers, and the solid cores are chemically etched to form an ordered array of pores comprised of lead silicate glass. Hydrogen firing is then used to activate the channel walls for electron multiplication. One drawback of this process is that the electrical conductivity and the secondary electron emission properties of the MCP cannot be adjusted independently, because both of these properties are imparted during the hydrogen-firing step. An additional drawback is that the aspect ratio (pore length to pore diameter) of the pores is limited to  $\sim 100$  because the chemical etchant cannot diffuse into longer pores without dissolving the lead silicate cladding. This limits the gain that can be achieved.

An alternative strategy for MCP manufacturing is to functionalize glass capillary array wafers using ALD to impart the necessary conductivity and secondary emission properties.<sup>15,16</sup> This process is illustrated in Figure 9 for a 1.2 mm thick, 33 mm diameter disc with a pore size of 20 μm (Figure 9A-C). Figure 9D shows the 33 mm capillary array disc after ALD of a 94 nm Mo:Al<sub>2</sub>O<sub>3</sub> composite film to serve as a resistive layer, and an 8 nm ALD MgO to enhance the secondary electron emission

of the surface.<sup>17</sup> After the ALD, nickel-chromium electrodes were evaporated onto both sides of the disc to generate the MCP shown in **Figure 9E**. The ability to tune the resistivity of the ALD composite film allows MCPs to be manufactured over a broad resistance range so as to be optimized for specific detector applications.

The fully functionalized MCP shown in **Figure 9E** was evaluated in a high vacuum system equipped with a calibrated electron source and a phosphor screen to determine its gain and spatial uniformity. As shown by the photograph of the phosphor screen image in **Figure 9F**, the MCP response is uniform across the surface. Next, a pair of identical MCPs was assembled in series and installed in front of a position-sensitive anode. Using a photon counting mode and weak ultraviolet illumination, a two-dimensional map of the spatial distribution of the gain was generated (**Figure 9G**), again showing the uniformity of the gain across the MCP. Studies on the performance of ALD MCPs have also been reported.<sup>14</sup>




**Figure 9.** MCP fabrication using ALD tunable resistance coatings. **A-C)** Increasing resolution of 33 mm diameter borosilicate glass capillary array substrate comprised of hexagonal bundles of 20  $\mu\text{m}$  pores. **D)** Glass substrate after 94 nm ALD Mo:Al<sub>2</sub>O<sub>3</sub> composite resistive layer and eight nm ALD MgO emissive layer. **E)** Completed MCP with nickel-chromium electrodes. **F)** Phosphor screen illuminated with electrons amplified by ALD MCP. **G)** Gain map for pair of ALD MCPs.

One of the benefits of using ALD for MCP manufacturing is that it facilitates scaling to larger detector sizes. The dark circular disk in **Figure 10** is the 33 mm MCP described above. To the left of this small MCP is an 8 x 8" glass capillary array substrate in a metal transport case. To the right is an identical 8 x 8" plate after ALD functionalization and electrode deposition to produce a large-area MCP. This is the largest MCP ever produced. Work is underway in the Large Area Picosecond Photodetector Project led by Argonne National Laboratory to implement these MCPs for applications in high-energy physics.<sup>14,15</sup>

### Charge Drain Coatings


Electrostatic charging is a challenging problem in electron-optical devices such as mass spectrometers, particle detectors, and micro-electro mechanical systems (MEMS). When stray electrons strike the surfaces of insulating components in a vacuum environment, the electrons become trapped, creating an electrostatic charge. With no means of draining the accumulated




**Figure 10.** Uncoated 8 x 8" glass capillary substrate in metal transport case (left) next to identical, 8x8" glass substrate after ALD Mo:Al<sub>2</sub>O<sub>3</sub> composite coating. Ruler and 33 mm MCP are shown for scale.

electrostatic charge, which can generate potentials of several kilovolts, the local electric fields become distorted causing loss of focus, resolution, and timing. Simple conducting or semiconducting films are unsuitable for draining this charge because the high leakage current can exceed the power supply or damage the device through resistive heating. Instead, a thin film with resistivity in the mesoscale (between those of a semiconductor and an insulator, 10<sup>2</sup>-10<sup>7</sup> Ohm cm), is needed. In electron-optical MEMS, the challenge is aggravated by the much closer proximity of the electrons to the surfaces, and the films are subjected to very high electrical fields of up to 25 MV/m. Furthermore, the film must coat all surface orientations while maintaining nanoscale uniformity to ensure equivalent performance across all components.

The ALD tunable resistance nanocomposite films are an excellent candidate for charge drain coatings. We have demonstrated their benefit in the reflective electron beam lithography (REBL) Nanowriter — a high throughput, direct-writing, maskless lithography tool under development at KLA—Tencor.<sup>18</sup> A schematic of the Nanowriter is shown in **Figure 11A**. The heart of the patterning system is the digital pattern generator (DPG) chip (**Figure 11B**). The DPG is an integrated circuit chip with an array of small, independently controllable electron mirrors ("lenslets") producing over 1 million "beamlets". The lenslet electrode MEMS structure is fabricated on top of complementary metal oxide semiconductor (CMOS) circuitry that provides voltage to four electrodes per lenslet to establish the electric fields required to produce an electrostatic lens for each pixel. The pixels and lenslets of the array are located on a 1.6  $\mu\text{m}$  pitch, and the lenslet structure is ~four  $\mu\text{m}$  in height. The DPG chip operates under a constant illumination of electrons, and this can lead to electrostatic charging on the insulating oxide surface along the wall of the lenslet (**Figure 11C**) that reduces the image quality (**Figure 12G**). To overcome this problem, a KLA-Tencor DPG chip was coated using an ALD Mo:Al<sub>2</sub>O<sub>3</sub> composite coating to drain static charge (**Figure 11D**). **Figures 12A-F** show images of electron beam test patterns obtained using this DPG. These images indicate virtually defect-free DPG patterns.<sup>19</sup>



**Figure 11.** A) Schematic of KLA-Tencor REBL Nanowriter. B) REBL digital pattern generator (DPG). C) Cleaved cross-section of DPG MEMS lenslet array structure. D) Schematic of ALD charge drain coating covering DPG lenslet array.



**Figure 12.** Defect-free test pattern images from DPG coated using ALD Mo:Al<sub>2</sub>O<sub>3</sub> composite charge drain coating. A) Square pattern, B) chevron pattern, C) stitching pattern, D) stig test pattern, E) KLA-Tencor company logo, and F) contact hole pattern. G) Poor quality pattern from DPG without ALD Mo:Al<sub>2</sub>O<sub>3</sub> charge drain coating.

### Solar Selective Coatings

Concentrated solar thermal power combines both energy harvesting and energy storage in a single plant, making it very appealing. Solar “power tower” facilities can generate hundreds of megawatts using a vast array of mirrors to reflect solar radiation onto a central receiver, where the energy is absorbed and used to heat a fluid that can be stored, or used to drive a turbine for electric power (Figure 13A). The efficiency can be increased using a solar-selective coating that absorbs strongly in the visible solar spectrum, while simultaneously reflecting infrared (IR) wavelengths. A high IR reflectivity ensures that the receiver will not radiate IR energy at high operating temperatures, allowing more of the energy to be captured and used. In addition, the coating must endure high temperatures without losing its selective absorption properties.


As described above, the optical absorption of the ALD W:Al<sub>2</sub>O<sub>3</sub> composite films can be tuned by adjusting the W content. This coating can be deposited on a metal substrate to absorb visible light while reflecting IR wavelengths. We prepared a 100 nm ALD W:Al<sub>2</sub>O<sub>3</sub> composite film using 33% W cycle percentage on a Hayes 230 nickel superalloy coupon that had been pre-coated

with a high surface area mesoporous film. The resulting material is opaque in the visible, but highly reflective in the IR. Figure 13B shows the substrate after 12 hours in air at 800 °C, and 50 thermal shock cycles between 300–800 °C in air. Most conventional materials would oxidize during this treatment and become transparent, but our ALD nanocomposite coating is extremely stable at high temperatures in air, and remains black. Moreover, our coating resists cracking and delamination during repeated thermal cycling that simulates the day-night temperature changes. These ALD materials are very promising for use as selective solar absorbing coatings in next-generation concentrated solar thermal power facilities.

### Conclusions

We have developed nanocomposite coatings built from conducting, metallic nanoparticles embedded in an amorphous dielectric matrix. These films are comprised of ALD M:D films where M= W, Mo, Ta, Nb or Co, and D= Al<sub>2</sub>O<sub>3</sub>, ZrO<sub>2</sub>, HfO<sub>2</sub>, MgO, or TiO<sub>2</sub>. By varying the ratio of ALD cycles for metal and Al<sub>2</sub>O<sub>3</sub> layers, we can tune precisely the resistance of these





**Figure 13.** A) Aerial photograph of Ivanpah concentrated solar power facility showing array of mirrors focussing sunlight on central receiving tower. B) 100nm ALD W:Al<sub>2</sub>O<sub>3</sub> composite films prepared using 33% W cycle percentage on Hayes 230 nickel superalloy coupon that had been pre-coated with high surface area mesoporous film after 12 hours in air at 800 °C, and 50 thermal shock cycles between 300–800 °C in air.

coatings over a very broad range from  $10^5$ – $10^{12}$  Ohm cm, and the bandgap and absorption coefficient across the visible spectrum. We employed *in situ* QCM and FTIR absorption spectroscopy to understand the ALD growth mechanism for the Mo:Al<sub>2</sub>O<sub>3</sub> nanocomposite films, showing that TMA serves as the reducing agent to form the conducting Mo. Cross-sectional TEM revealed the film structure to be metallic nanoparticles (1–2 nm) embedded in an amorphous matrix. We utilized these nanocomposite coatings to functionalize capillary glass array plates to fabricate large-area MCPs suitable for application in large-area photodetectors. In addition, we have applied these films as charge drain coatings in MEMS devices for a prototype electron beam lithography tool, and as selective solar absorbing coatings. These practical applications promise a bright future for the use of ALD for creation of nanocomposite coatings for specific commercial applications.

## Acknowledgements

This work was supported by the U.S. Department of Energy, Office of Science, Office of Basic Energy Sciences and Office of High Energy Physics under contract DE-AC02-06CH11357 as part of the Large Area Picosecond Photodetector (LAPPD) project. The work at KLA-Tencor was partly sponsored by Defense Advanced Research Projects Agency under contract number HR0011-07-9-0007. The views, opinions, and/or findings contained in this article/presentation are those of the author/presenter and should not be interpreted as representing the official views or policies, either expressed or implied, of the Defense Advanced Research Projects Agency or the Department of Defense.














## References








- (1) George, S. M. *Chem. Rev.* **2010**, *110*(1), 111–113.
- (2) Elam, J. Coatings on High Aspect Ratio Structures in *Atomic Layer Deposition of Nanostructured Materials*. Pinna, N. and Knez, M., Eds. Weinheim: Wiley-VCH. **2012**, 227.
- (3) Elam, J. W.; George, S. M. *Chem. Mater.* **2003**, *15*, 1020.
- (4) Elam, J. W.; Routkevitch, D.; George, S. M. *J. Electrochem. Soc.* **2003**, *150*, 6, G339.
- (5) Mane, A. U.; Elam, J. W. *Chem. Vap. Deposition* **2013**, *19*, 186.
- (6) Elam, J. W.; Groner, M. D.; George, S. M. *Rev. Sci. Instrum.* **2002**, *73*(8), 2981.
- (7) Brodie, A.; De Cecco, P.; Bevis, C.; Maldonado, J. R.; Bhatia, R.; Deguns, E.; Sundaram, G. Fabrication of Coatings with Targeted Tunable Electrical Properties via ALD: Al<sub>2</sub>O<sub>3</sub>/ZnO and Nb<sub>2</sub>O<sub>5</sub>/Ta<sub>2</sub>O<sub>5</sub> in *Atomic Layer Deposition Applications 6*. Elam, J. W.; DeGendt, S.; VanDerStraten, O.; Delabie, A.; Londergan, A.; Bent, S. F.; Roozeboom, F.; Eds. **2010**, *33*, 101.
- (8) Ferguson, J. D.; Weimer, A. W.; George, S. M. *Appl. Surf. Sci.* **2000**, *180*, 162–163.
- (9) Seghete, D.; Rayner, Jr. G. B.; Cavanagh, A. S.; Anderson, V. R.; George, S. M. *Chem. Mater.* **2011**, *23*(7) 1668.
- (10) Christensen, S. T.; Elam, J. W. *Chem. Mater.* **2010**, *22*(8), 2517.
- (11) Du, X.; Du, Y.; George, S. M. *J. Phys. Chem. A*, **2008**, *112*(39) 9211.
- (12) Klug, J. A.; Proslie, T.; Elam, J. W.; Cook, R. E.; Hiller, J. M.; Claus, H.; Becker, N. G.; Pellin, M. J. *J. Phys. Chem. C*, **2011**, *115*(50) 25063.
- (13) Fabreguette, F. H.; Sechrist, Z. A.; Elam, J. W.; George, S. M. *Thin Solid Films* **2005**, *488*(1–2), 103.
- (14) Siegmund, O. H. W.; McPhate, J.; Tremsin, A. S.; Jelinsky, S. R.; Frisch, H. J.; Elam, J. W.; Mane, A. U.; Wagner, R. G.; Minot, M. J.; Renaud, J. M.; Detarando, M. A. *Phys. Procedia* **2012**, *37*, 803.
- (15) Mane, A. U.; Peng, Q.; Elam, J. W.; Bennis, D. C.; Craven, C. A.; Detarando, M. A.; Escolás, J. R.; Frisch, H. J.; Jokela, S. J.; McPhate, J.; Minot, M. J.; Siegmund, O. H. W.; Renaud, J. M.; Wagner, R. G.; Wetstein, M. J. An atomic layer deposition method to fabricate economical and robust large area microchannel plates for photodetectors, in *Proceedings of the 2nd International Conference on Technology and Instrumentation in Particle Physics*. Liu, T.; Ed. **2012**, *37*, 722.
- (16) D. R. Beaulieu, D. Gorelikov, H. Klotzsch, P. de Rouffignac, K. Saadatmand, K. Stenton, N. Sullivan, and A. S. Tremsin, *Nucl. Instrum. Methods Phys. Res. A* **2011**, *633*, S59.
- (17) Jokela, S. J.; Veryovkin, I. V.; Zinovev, A. V.; Elam, J. W.; Mane, A. U.; Peng, Q.; Insepov, Z.; et al. Secondary electron yield of emissive materials for large-area micro-channel plate detectors: surface composition and film thickness dependencies in *Proceedings of the 2nd International Conference on Technology and Instrumentation in Particle Physics*. Liu, T.; Ed. **2012**, *37*, 740.
- (18) Petric, P.; Bevis, C.; McCord, M.; Carroll, A.; Brodie, A.; Ummethala, U.; Grella, L.; Cheung, A.; Freed, R. *J. Vac. Sci. Technol. B* **2010**, *28*(6), C6C6.
- (19) Tong, W. M.; Brodie, A. D.; Mane, A. U.; Fuge, S.; Kidwingira, F.; McCord, M. A.; Bevis, C. F.; Elam, J. W. *Appl. Phys. Lett.* **2013**, *102*(25), 252901.



## Precursors Packaged for Deposition Systems

For a complete list of available materials, visit [SigmaAldrich.com/ald](http://SigmaAldrich.com/ald).

Name	Acronym	Structure	Form	Cat. No.
(3-Aminopropyl)triethoxysilane	APTES; APTS		liquid	706493-20ML
Bis(ethylcyclopentadienyl)cobalt(II)	(EtCp) <sub>2</sub> Co		liquid	753076-10G
Bis(methyl-η <sup>5</sup> -cyclopentadienyl)methoxymethylzirconium	ZrD-CO <sub>4</sub> ; ZRCMMM		liquid	725471-10G
Diethylzinc	DEZ; DEZn; Et <sub>2</sub> Zn		liquid	668729-25G
Niobium(V) ethoxide	NbOEt	$\text{Nb}(\text{OCH}_2\text{CH}_3)_5$	liquid	760412-10G
Silicon tetrachloride	STC	$\text{SiCl}_4$	liquid	688509-25ML
Tetrakis(diethylamido)titanium(IV)	TDEAT		liquid	725536-10G
Tetrakis(dimethylamido)hafnium(IV)	TDMAH		low-melting solid	666610-25G
Tetrakis(dimethylamido)titanium(IV)	TDMAT		liquid	669008-25G
Tetrakis(dimethylamido)zirconium(IV)	TDMAZ		solid	669016-25G
Tetrakis(ethylmethylamido)hafnium(IV)	TEMAH		liquid	725544-10G
Tetrakis(ethylmethylamido)zirconium(IV)	TEMAZ		liquid	725528-10G
Titanium(IV) isopropoxide	TTIP		liquid	687502-25G
Titanium tetrachloride	TTC	$\text{TiCl}_4$	liquid	697079-25G
Triethylgallium	TEG; TEGa; Et <sub>3</sub> Ga		liquid	730726-10G
Trimethylaluminum	TMA		liquid	663301-25G

Name	Acronym	Structure	Form	Cat. No.
Trimethylgallium	TMG; TMGa; Me <sub>3</sub> Ga		liquid	730734-10G
Trimethyl(methylcyclopentadienyl)platinum(IV)	MeCpPtMe <sub>3</sub>		low-melting solid	697540-10G
Tris[N,N-bis(trimethylsilyl)amide]yttrium	YDTMSA		powder	702021-10G
Tris(diethylamido)(tert-butylimido)niobium(V)	-		liquid	751774-10G
Tris(diethylamido)(tert-butylimido)tantalum(V)	TBTDET		liquid	668990-10G
Tris(dimethylamino)silane	TDMAS		liquid	759562-25G
Tungsten hexacarbonyl	-	W(CO) <sub>6</sub>	solid	755737-25G
Water	-	H <sub>2</sub> O	liquid	697125-25ML
Zirconium(IV) tert-butoxide	ZTB		liquid	759554-25G

## Sputtering Targets

For a complete list of available materials, visit [SigmaAldrich.com/sputtering-target](https://www.sigmaaldrich.com/sputtering-target).

Name	Structure	Purity (%) Or Grade	Cat. No.
Aluminum	Al	99.9995% trace metals basis, sputtering target	749036-1EA
Aluminum zinc oxide	Al <sub>2</sub> O <sub>3</sub> / ZnO	99.99% trace metals basis, sputtering target	752665-1EA
Chromium	Cr	99.95% trace metals basis, sputtering target	749052-1EA
Copper	Cu	99.95% trace metals basis, sputtering target	767476-1EA
Indium tin oxide	In <sub>2</sub> O <sub>3</sub> /SnO <sub>2</sub>	99.99% trace metals basis, sputtering target	752657-1EA
Nickel	Ni	99.95% trace metals basis, sputtering target	767484-1EA
Silicon	Si	99.999% trace metals basis, sputtering target	767492-1EA
Tantalum	Ta	99.95% trace metals basis, sputtering target	767514-1EA
Titanium	Ti	99.995% trace metals basis, sputtering target	749044-1EA
		99.995% trace metals basis, sputtering target	767506-1EA
Tungsten	W	99.95% trace metals basis, sputtering target	767522-1EA
Yttrium(III) oxide	Y <sub>2</sub> O <sub>3</sub>	99.99% trace metals basis, sputtering target	774022-1EA
Yttrium	Y	99.9% trace metals basis, sputtering target	773972-1EA
Zinc oxide	ZnO	99.99% trace metals basis, sputtering target	752681-1EA
Zinc	Zn	99.995% trace metals basis, sputtering target	749060-1EA
Zirconium(IV) oxide	ZrO <sub>2</sub>	99.95% trace metals basis (excludes 2% HfO <sub>2</sub> ), sputtering target	774030-1EA
Zirconium(IV) oxide-yttria stabilized	ZrO <sub>2</sub> / Y <sub>2</sub> O <sub>3</sub>	99.9% trace metals basis, sputtering target	774049-1EA
Zirconium yttrium alloy	Zr <sub>0.85</sub> Y <sub>0.15</sub>	99.9% trace metals basis (excluding ≤1% Hf), sputtering target	774057-1EA

# Silicon Nitride Atomic Layer Deposition: A Brief Review of Precursor Chemistry



Antonio T. Lucero, Jiyoung Kim\*

Department of Materials Science and Engineering,  
The University of Texas at Dallas, Richardson, TX 75080, USA  
\*Email: jiyoung.kim@utdallas.edu

## Introduction

Silicon nitride ( $\text{SiN}_x$ ) is a critical material for semiconductor devices, increasingly used in high-performance logic and memory. Modern, scaled devices require robust  $\text{SiN}$  films deposited at low temperature ( $<400$  °C) for use as gate sidewall spacers and in self-aligned quadruple patterning.<sup>1</sup> Traditional  $\text{SiN}_x$  deposition techniques, including chemical vapor deposition (CVD) and plasma-enhanced chemical vapor deposition (PECVD), are now giving way to atomic layer deposition (ALD). ALD allows more control over the thickness of deposition, work at relatively low temperatures, and conforms over high-aspect ratio structures.<sup>2</sup> ALD can be divided into two classes, thermal ALD and plasma-enhanced ALD (PEALD). Both methods have some advantages for  $\text{SiN}_x$  deposition. Thermal ALD allows for conformal deposition over high aspect ratio (HAR) structures ( $>5000:1$ ), while PEALD can be used at much lower temperatures with lower HAR conformality. Advances in precursor chemistry and nitrogen sources have enabled the tailoring of material properties like the wet etch rate and growth rate to meet research and industry requirements. There are currently three main silicon precursor classes: chlorosilanes, organosilanes, and heterosilanes. Chlorosilanes are silicon precursors where Si-Cl bonding is predominant. Organosilanes are silicon precursors containing organic ligands, although currently this class is limited to aminosilanes in practice. The last group, heterosilanes, includes all other precursors.

## Chlorosilanes

Chlorosilanes are an historically important class of silicon precursors that helped to build the semiconductor industry by enabling the production of ultra-high purity silicon. This class includes any silicon precursor containing at least one chlorine-silicon bond. The first  $\text{SiN}_x$  was grown by thermal ALD in 1997, when Morishita<sup>3</sup> deposited  $\text{SiN}_x$  using hexachlorodisilane (HCDS,  $\text{Si}_2\text{Cl}_6$ , **Cat. No. 205184**) and hydrazine ( $\text{N}_2\text{H}_4$ , **Cat. No. 215155**) at temperatures ranging from 525–650 °C. While hydrazine has since been replaced by more convenient nitrogen sources, HCDS has remained an important precursor for  $\text{SiN}_x$  ALD. Later reports<sup>4</sup> of ALD using HCDS and ammonia demonstrate successful deposition of  $\text{SiN}_x$  at temperatures as

low as 515–557 °C. In addition, tetrachlorosilane,<sup>5</sup> dichlorosilane (DCS),<sup>6</sup> and octachlorotrisilane<sup>7</sup> have all been used successfully. Since the deposition temperature is fairly high, the physical properties such as density and wet etch rate (WER) in hydrofluoric acid are good. Growth per cycle (GPC) varies, but is typically greater than 1 Å/cycle. A disadvantage of chlorosilanes for thermal ALD of  $\text{SiN}_x$  is the large precursor exposure ( $10^7$ – $10^{10}$  L) required to achieve saturation. Of note, only chlorosilane precursors have been used for thermal ALD of  $\text{SiN}_x$ , since they are the only precursors stable enough for use above 400 °C, the temperature at which ammonia is activated. With the exception of hydrazine, which will be discussed in a later section, no other nitrogen sources are available for thermal ALD. This limitation has prevented the widespread industry adoption of thermal ALD for  $\text{SiN}_x$  growth.

In order to grow films at low ( $<400$  °C) temperatures, research has focused on the use of plasma to aid deposition.<sup>8</sup> Microwave plasma, inductively coupled plasma (ICP), and capacitively coupled plasma (CCP) are most commonly used in combination with sources of reactive nitrogen including ammonia, nitrogen, or nitrogen forming gas ( $\text{N}_2$ - $\text{H}_2$ ). DCS<sup>9</sup> and HCDS<sup>10</sup> have been used extensively with ammonia in PEALD of  $\text{SiN}_x$  at temperatures from 300–400 °C. Temperatures less than 300 °C can lead to excess chlorine contamination due to the formation of  $\text{NH}_4\text{Cl}$ . Ovanesyan et al.<sup>10</sup> reported conformal deposition of  $\text{SiN}_x$  on HAR structures using HCDS and  $\text{NH}_3$  plasma at 400 °C, with hydrogen in the form of -NH and chlorine ( $<1\%$ ) as the primary impurities. Conformal deposition when using  $\text{NH}_3$  plasma is a strong advantage of chlorosilane precursors. Unfortunately, the WER for  $\text{SiN}_x$  deposited with chlorosilane precursors is reported to be high, and film density is low due to hydrogen incorporation. Recently, a new chlorosilane precursor, pentachlorodisilane (PCDS) has been reported<sup>11</sup> that while similar to HCDS, results in a 20% higher GPC (0.78 vs. 1.02 Å/cycle) with similar or better physical properties. Substituting a chlorine atom with a hydrogen appears to lower the steric hindrance of the PCDS molecule and increase its polarity, leading to a precursor with higher reactivity. In addition, a precursor exposure of only  $4 \times 10^4$  L,

or 4–5 orders of magnitude lower than the exposure for thermal ALD processes and 1–2 orders lower than that of other PEALD provides results with a similar GPC. The unique hollow cathode plasma source used for films grown with both HCDS and PCDS results in exceptionally low oxygen contamination in these films.

### Organosilanes

The first organosilane used for SiN<sub>x</sub> ALD was tris(dimethylamino) silane (TDMAS, **Cat. Nos. 570133, 759562**) in 2008.<sup>12</sup> Using a remote ICP nitrogen-forming gas plasma to successfully deposit SiN<sub>x</sub>, though with carbon impurities of 5–10%. Provine et al.<sup>13</sup> improved on these results and were able to grow SiN<sub>x</sub> with high film density (2.4 g/cm<sup>3</sup>) and low WER (3 nm/min in 100:1 HF) at a temperature of 350 °C. Performing a hydrogen plasma post-anneal reduced WER to less than 1 nm/min.

Bis(tert-butylamino)silane (BTBAS) is another aminosilane frequently used for SiN<sub>x</sub> deposition. Knoops et al. deposited high quality SiN<sub>x</sub> with BTBAS and N<sub>2</sub> plasma.<sup>14</sup> Film density was very high at 2.8 g/cm<sup>3</sup>, and film wet etch rate was 0.2 nm/min for growth at 400 °C. Carbon contamination was less than 2%, but that increased to approximately 10% for films grown at 200 °C. Film properties were similar to those obtained from low-pressure chemical vapor deposition (LPCVD) grown SiN<sub>x</sub>, which is attributed to the high film density of the film.

All organosilane precursors using nitrogen plasma, causes the GPC to drop to nearly zero. When NH<sub>3</sub> plasma is used, -NH<sub>2</sub> is a common surface termination. Huang et al.<sup>15</sup> predicted a high activation energy barrier for BTBAS and other amine groups, and this prediction was experimentally confirmed by Provine for TDMAS and NH<sub>3</sub> plasma growth. Density functional theory simulations predicted that BTBAS would only adsorb at undercoordinated nitrogen or silicon sites, which NH<sub>3</sub> plasma cannot provide. Another issue with most organosilanes is the low GPC — typically less than 0.3 Å/cycle. Organosilanes are large molecules that provide only a single silicon atom, so steric hindrance likely plays a role, though the requirement for undercoordinated surface sites probably hinders growth as well. Finally, conformality over HAR structures is limited for organosilanes. Faraz et al. found at best 50% conformality for sidewalls (26 vs 13 nm) when using di(sec-butylamino)silane (DSBAS) and N<sub>2</sub> plasma.<sup>1</sup> These results are common, and are likely due to the soft saturation observed for N<sub>2</sub> plasma exposure. Typically sidewalls are exposed to lower plasma density during the deposition process compared to the top and bottom of the structure.

### Heterosilanes

The last group of precursors encompasses all non-organic and non-halide precursors. Silica is simplest of the precursors and has already been discussed. Well-known for use its use in PECVD, SiN<sub>x</sub> can be deposited using SiH<sub>4</sub> and N<sub>2</sub> plasma.<sup>16</sup> The disadvantage of this approach is the long (60 s) N<sub>2</sub> plasma time required for saturation. Likewise, S-H contamination is problematic and likely results in low film density and high WER. Trisilylamine (TSA) is another Si precursor that has been

investigated for SiN<sub>x</sub> deposition. Triyoso et al.<sup>17</sup> demonstrated the growth of SiN<sub>x</sub> with TSA and N<sub>2</sub>/H<sub>2</sub> plasma at 300 and 400 °C. GPC depended greatly on plasma conditions, and varied from 1.3 to 2.1 Å/cycle and WER in 100:1 HF was approximately 1 nm/min. When comparing PEALD SiN<sub>x</sub> with PECVD SiN<sub>x</sub>, they found TSA-based PEALD improved transistor performance with higher drive current, higher hole mobility, and improved I<sub>on</sub>/I<sub>off</sub> ratio. Jang et al.<sup>18</sup> deposited SiN<sub>x</sub> with TSA and NH<sub>3</sub> plasma, though with a lower GPC (0.6 Å/cycle vs. ~1.5 Å/cycle). This has important implications for conformal deposition over HAR structures, since NH<sub>3</sub> plasma offers improved step coverage. Finally, neopentasilane (NPS) and N<sub>2</sub> plasma were used to grow SiN<sub>x</sub> between 250 and 300 °C.<sup>19</sup> Film properties were found to be similar to that of TSA, though NPS has a slightly higher GPC (1.2 vs. 1.4 Å/cycle). The WER for SiN<sub>x</sub> grown with NPS was strongly plasma dependent, but optimized conditions resulted in a WER of between 2 and 3 nm/min. Both TSA and NPS are interesting precursors due to their high silicon percentage relative to molecular weight, and they offer higher GPC than organosilanes while maintaining a lower WER than typical chlorosilanes.

### Summary and Outlook

The need for high quality, conformal SiN<sub>x</sub> films grown at low temperature is increasing, and both academia and industry are working to develop advanced processes and precursors. Currently available precursors offer a range of advantages and disadvantages. Chlorosilanes provide enhanced GPC and conformality over HAR structures, but lack wet etch resistance and film density. Organosilanes enable the growth of SiN<sub>x</sub> films with extremely low WER, comparable to or lower than those achieved with LPCVD, but are hindered by low GPC and poor conformality. The heterosilanes, which include trisilylamine and neopentasilane, offer both good GPC, which aids in throughput, and low WER. Further investigation is needed to see if these precursors can deliver conformality over HAR structures.

While not a silicon precursor, it is worth discussing current work on nitrogen sources. Recent advances in hydrazine delivery technology have enabled the use of ultra-high purity hydrazine in thermal ALD. While toxicity is still a concern, the overall safety of the source has been improved.<sup>20</sup> Thermal ALD of TaN and WN with hydrazine as the nitrogen source has been demonstrated. Likewise, low resistivity of TiN deposited at 275–350 °C has been reported. Few reports exist regarding the deposition of SiN<sub>x</sub> with hydrazine, but the growth of SiN<sub>x</sub> passivation layers using HCDS and hydrazine at 285 °C has been successful.<sup>21</sup> Chlorine contamination is a problem; these issues can likely be associated with the extremely low deposition temperature since similar issues occur for PEALD below 300 °C. Current work using ultra-high purity hydrazine as a source for SiN<sub>x</sub> deposition is promising, and good SiN<sub>x</sub> properties have been obtained in depositions occurring between 350 and 400 °C, which should enable the development of low temperature thermal ALD of SiN<sub>x</sub>. Together, improvements in nitrogen sources and silicon precursors portend a bright future for atomic layer deposition of SiN<sub>x</sub>.












## References

- (1) Faraz, T.; van Druenen, M.; Knoops, H. C.; Mallikarjunan, A.; Buchanan, I.; Hausmann, D. M.; Henri, J.; Kessels, W. M. *ACS Appl. Mater. Interfaces* **2017**, *9* (2), 1858–1869.
- (2) Meng, X.; Byun, Y. C.; Kim, H. S.; Lee, J. S.; Lucero, A. T.; Cheng, L.; Kim, J. *Materials (Basel)* **2016**, *9* (12).
- (3) Morishita, S.; Sugahara, S.; Matsumura, M. *Appl. Surf. Sci.* **1997**, *112*, 198–204.
- (4) Park, K.; Yun, W.-D.; Choi, B.-J.; Kim, H.-D.; Lee, W.-J.; Rha, S.-K.; Park, C. O. *Thin Solid Films* **2009**, *517* (14), 3975–3978.
- (5) Nakajima, A.; Yoshimoto, T.; Kidera, T.; Obata, K.; Yokoyama, S.; Sunami, H.; Hirose, M. *Appl. Phys. Lett.* **2000**, *77* (18), 2855–2857.
- (6) Hansch, W.; Nakajima, A.; Yokoyama, S. *Appl. Phys. Lett.* **1999**, *75* (11), 1535–1537.
- (7) Riedel, S.; Sundqvist, J.; Gumprecht, T. *Thin Solid Films* **2015**, *577*, 114–118.
- (8) Profijt, H. B.; Potts, S. E.; van de Sanden, M. C. M.; Kessels, W. M. M. *J. Vac. Sci. Technol. A* **2011**, *29* (5), 050801.
- (9) Nagata, K.; Nagasaka, M.; Yamaguchi, T.; Ogura, A.; Oji, H.; Son, J.-Y.; Hirotsawa, I.; Watanabe, Y.; Hirota, Y. *ECS Transactions* **2013**, *53* (3), 51–56.
- (10) Ovanesyan, R. A.; Hausmann, D. M.; Agarwal, S. *ACS Appl. Mater. Interfaces* **2015**, *7* (20), 10806–10813.
- (11) Meng, X.; Kim, H. S.; Lucero, A. T.; Hwang, S. M.; Lee, J. S.; Byun, Y.-C.; Kim, J.; Hwang, B. K.; Zhou, X.; Young, J.; Telgenhoff, M. *ACS Appl. Mater. Interfaces* **2018**, *10* (16), 14116–14123.
- (12) F. Q. and H. C., in *Silicon Nitride and Silicon Oxide Thin Films by Plasma ALD*, Bruges, Belgium, **2008** (Proceedings of the 8th International Conference on Atomic Layer Deposition).
- (13) Provine, J.; Schindler, P.; Kim, Y.; Walch, S. P.; Kim, H. J.; Kim, K.-H.; Prinz, F. B. *AIP Adv.* **2016**, *6* (6), 065012.
- (14) Knoops, H. C.; Braeken, E. M.; de Peuter, K.; Potts, S. E.; Haukka, S.; Pore, V.; Kessels, W. M. *ACS Appl. Mater. Interfaces* **2015**, *7* (35), 19857–62.
- (15) Huang, L.; Han, B.; Han, B.; Derecskei-Kovacs, A.; Xiao, M.; Lei, X.; O'Neill, M. L.; Pearlstein, R. M.; Chandra, H.; Cheng, H. *Phys. Chem. Chem. Phys.* **2014**, *16* (34), 18501–18512.
- (16) King, S. W. *J. Vac. Sci. Technol. A* **2011**, *29* (4), 041501.
- (17) Triyoso, D. H.; Hempel, K.; Ohsiek, S.; Jaschke, V.; Shu, J.; Mutas, S.; Dittmar, K.; Schaeffer, J.; Utess, D.; Lenski, M. *ECS J. Solid State Sci. Technol.* **2013**, *2* (11), N222–N227.
- (18) Jang, W.; Jeon, H.; Song, H.; Kim, H.; Park, J.; Kim, H.; Jeon, H. *Phys. Status Solidi A* **2015**, *212* (12), 2785–2790.
- (19) Weeks, S.; Nowling, G.; Fuchigami, N.; Bowes, M.; Littau, K. *J. Vac. Sci. Technol. A* **2016**, *34* (1), 01A140.
- (20) Alvarez, D.; Spiegelman, J.; Holmes, R.; Andachi, K.; Raynor, M.; Shimizu, H. *ECS Transactions* **2017**, *77* (5), 219–225.
- (21) Edmonds, M.; Sardashti, K.; Wolf, S.; Chagarov, E.; Clemons, M.; Kent, T.; Park, J. H.; Tang, K.; McIntyre, P. C.; Yoshida, N.; Dong, L.; Holmes, R.; Alvarez, D.; Kummel, A. C. *J. Chem. Phys.* **2017**, *146* (5), 052820.

## Silicon ALD Precursors

For a complete list of available materials, visit [SigmaAldrich.com/vapordeposition](http://SigmaAldrich.com/vapordeposition).












Name	Structure	Purity (%) Or Grade	Form	Cat. No.
(3-Aminopropyl)triethoxysilane		≥98.0%	liquid	<b>741442-100ML</b> <b>741442-500ML</b>
<i>N</i> -sec-Butyl(trimethylsilyl)amine		98%	liquid	<b>654671-5G</b>
Chloropentamethylidisilane		97%	liquid	<b>490407-5G</b>
1,2-Dichlorotetramethyldisilane		95%	liquid	<b>415456-5ML</b>
1,3-Diethyl-1,1,3,3-tetramethyldisilazane		>98.0%	liquid	<b>646261-1G</b> <b>646261-5G</b>
Dodecamethylcyclohexasilane		≥97%, GC	powder or crystals	<b>437492-1G</b>
Hexamethyldisilane		98%	liquid	<b>217069-5G</b> <b>217069-10G</b> <b>217069-50G</b>
Hexamethyldisilazane		99.9%	liquid	<b>379212-25ML</b> <b>379212-100ML</b>
		≥99%	liquid	<b>440191-100ML</b> <b>440191-1L</b>
Methylsilane		≥99.9%	liquid	<b>462993-10G</b> <b>462993-20G</b>

Name	Structure	Purity (%) Or Grade	Form	Cat. No.
2,4,6,8,10-Pentamethylcyclopentasiloxane		96%	liquid	<b>517801-25ML</b>
Pentamethyldisilane		97%	liquid	<b>557641-5ML</b>
Silicon tetrabromide	SiBr <sub>4</sub>	99.995% trace metals basis	liquid	<b>494100-5ML</b> <b>494100-25ML</b>
Silicon tetrachloride	SiCl <sub>4</sub>	99%	liquid	<b>333468-10G</b>
		99.998% trace metals basis, packaged for use in deposition systems	liquid	<b>289388-100ML</b> <b>289388-800ML</b>
		99%	liquid	<b>688509-25ML</b> <b>215120-100G</b> <b>215120-1KG</b>
Tetraethylsilane		99%	liquid	<b>510874-5ML</b>
2,4,6,8-Tetramethylcyclotetrasiloxane		≥99.999% trace metals basis ≥98.5%	liquid	<b>512990-25ML</b> <b>512990-100ML</b>
1,1,2,2-Tetramethyldisilane		98%	liquid	<b>445126-1G</b>
Tetramethylsilane		≥99.99% trace metals basis	liquid	<b>523771-100ML</b>
Tris( <i>tert</i> -butoxy)silanol		packaged for use in deposition systems	solid	<b>697281-25G</b>
Tris( <i>tert</i> -pentoxy)silanol		packaged for use in deposition systems	liquid	<b>697303-25G</b>

## Silicon Solution Deposition Precursors

For a complete list of available materials, visit [SigmaAldrich.com/solutiondeposition](https://www.sigmaaldrich.com/solutiondeposition).

Name	Structure	Purity (%) Or Grade	Cat. No.
Allyltrimethoxysilane		≥98%	<b>679267-50G</b>
Butyltrichlorosilane		≥98%	<b>679224-50G</b>
[3-(Diethylamino)propyl]trimethoxysilane		≥98%	<b>679356-50G</b>

Name	Structure	Purity (%) Or Grade	Cat. No.
Dimethoxydimethylsilane		≥99.5% 99.999% metals basis	556688-25ML
Isobutyl(trimethoxy)silane		≥98%	679364-50G
Methyltrichlorosilane		≥98%, GC ≥99.99% (as metals)	679208-50G
<i>n</i> -Propyltriethoxysilane		≥98%	679321-50G
Triethoxy(octyl)silane		99.99% trace metals basis 98%	679305-50G
Triethoxyphenylsilane		≥98%	679291-50G
Triethoxyvinylsilane		≥98%	679275-50G
Trimethoxymethylsilane		≥98%	679232-50G
Trimethoxyphenylsilane		98%	679313-50G
Tris( <i>tert</i> -butoxy)silanol		99.999%	553468-5G 553468-25G
Tris( <i>tert</i> -pentoxy)silanol		≥99.99%	553441-5G 553441-25G

# Group 11 Thin Films by Atomic Layer Deposition




Sean Barry

Department of Chemistry, Carleton University, Ottawa,  
Ontario K1S 5B6, Canada  
Email: sean\_barry@carleton.ca

## Introduction

Copper metal deposition processes are an essential tool for depositing interconnects used in microelectronic applications, giving group 11 (coinage metals: Copper, Silver, and Gold) an important place in atomic layer deposition (ALD) process development. A significant amount of development has been invested in the design and improvement of both processes and precursors of  $\text{Cu}^+$  and  $\text{Cu}^{2+}$  oxidation states over the last 20 years (see Figure 1).



**Figure 1.** The cumulative number of coinage metal ALD processes developed over the past 20 years.

The timeline for the development of silver and gold metal deposition precursors roughly follows their standard potentials (Table 1). This is probably due to ease of synthesis and isolation, since the reactivity of metal deposition precursors is loosely governed by their standard potential. A key element of precursor development is testing thermal stability, so processes involving  $\text{Ag}^+$  and  $\text{Au}^{3+}$  were quickly developed, but many high potential cations have still not been successfully incorporated into precursors.

**Table 1.** The standard potentials (against a hydrogen electrode) for the common cations of the coinage metals.


Redox Half Reaction	Standard Potential, $E^0$ (V)
$\text{Ag}^{2+} + 2 e^- \rightarrow \text{Ag}^0$	1.98
$\text{Au}^+ + e^- \rightarrow \text{Au}^0$	1.69
$\text{Au}^{3+} + 3 e^- \rightarrow \text{Au}^0$	1.40
$\text{Ag}^+ + e^- \rightarrow \text{Ag}^0$	0.80
$\text{Cu}^+ + e^- \rightarrow \text{Cu}^0$	0.52
$\text{Cu}^{2+} + 2 e^- \rightarrow \text{Cu}^0$	0.34

In this article, we will discuss coinage metal deposition processes in order to provide a sense of the most critical precursors, reducing agents, and processes. This review is not comprehensive, but rather a brief perspective on the current state of the art. For a more in-depth discussion on precursor and process tabulation see the recent review by Puurunen.

## Copper


Copper metal has been very well studied and a tremendous number of copper-based precursors have been developed, including halides, amidinates, guanidinates,  $\beta$ -diketonates, and aminoalkoxides (Figure 2).<sup>1</sup> In general, thermal copper depositions can occur over a wide temperature range (120–500 °C), while plasma and catalyzed depositions occur as low as room temperature, with the energy necessary for the process being supplied by the plasma, or circumvented by the catalyst. Interestingly, many depositions of copper by ALD start to exhibit non-self-limited, continuous growth (i.e., chemical vapor deposition, CVD) at 200–230 °C. Precursor examples (Figure 2) include  $\text{Cu}^+$  centers but are dominated by  $\text{Cu}^{2+}$ .





**Figure 2.** A representative array of currently available precursors for copper metal deposition.


The  $\text{Cu}^{2+}$  dication is reduced relatively easily to copper metal measured against the oxidation of  $\text{H}_2$ , but  $\text{Cu}^{2+}$  is the most difficult to reduce of any of the common oxidation states of any coinage metal. This might account for the tendency of copper precursors to have chemisorbed moieties containing  $\text{Cu}^+$ , even if the precursor contained  $\text{Cu}^{2+}$  (see **Table 1**). A good example of this phenomenon can be seen with the copper(II) aminoalkoxides (**Figure 3**).



**Figure 3.** Thermal dehydrogenation of an aminoalkoxide at a surface copper.

The mechanism for the thermal dehydrogenation of an aminoalkoxide at a copper surface was first elucidated in 1993 and the model for chemisorption was published in 2012. The model described chemisorption of the  $\text{Cu}^+$  species, while experimental work completed in 1996 demonstrated the presence of both the aldehyde and the alcohol from the aminoalkoxide ligand. This thermal behavior is not self-limiting, yet Bis(dimethylamino-2-propoxy)copper(II) ( $\text{Cu}(\text{DMAP})_2$ ) has been shown to undergo ALD using a variety of reducing agents including diethyl zinc, formic acid/hydrazine, and

boranedimethylamine. The additional methyl group in the beta position prevents aldehyde formation and arrests thermal decomposition. The DMAP precursor is still able to undergo ALD growth, including by reaction with amino borane, which is essentially a source of hydrogen at the reaction temperature to produce copper metal at the surface (**Figure 4**).<sup>7</sup>



**Figure 4.** A copper DMAP surface moiety can be reduced to metallic copper and the parent ligand by dihydrogen generated from the thermal decomposition of borane dimethylamine.


An interesting note about  $\text{Cu}(\text{DMAP})_2$  is that it has covalently bonded oxygen. From a precursor design point of view, this bonding arrangement is generally avoided since any atom bonded directly to the metal center is at considerable risk of becoming incorporated in the metal film as an impurity. Obviously, in the case of this precursor, this is unfounded.

## Silver

Silver metal ALD processes are much scarcer due to silver's lower standard potential and thus lower reactivity. It can be difficult to design a silver-containing precursor that has the thermal (and photo) stability to withstand volatilization. The first silver metal deposited by ALD was in 2007, from a phosphine-supported silver(I) carboxylate. This process had a growth rate of 1.2 Å/cycle, but was only viable at 140 °C due to both volatility and a low thermal decomposition temperature of the precursor. Additionally, significant impurities were observed in the resulting film, including 10% oxygen contamination.


In 2014, trimethylphosphine-supported silver(I) hexafluoroacetoacetate was shown to undergo thermal ALD of silver metal using both formalin (aqueous formaldehyde stabilized with methanol) in a two step type process, and with trimethylaluminum (TMA, **Cat. No. 663301**) and water in three step process. In the first instance, the silver metal had a growth rate of 0.7 Å/cycle at 200 °C, although it showed somewhat slower growth at 170 °C. With TMA and water, the growth rate was significantly lower, about 0.2 Å/cycle at 110 °C. Interestingly, the growth rate could be improved with multiple, sequential doses of TMA/water, which appeared not to show saturation. The authors speculated that the TMA and water pulses were removing hexafluoroacetoacetate moieties from the surface, allowing better nucleation of the incoming precursor.

The precursor hexafluoroacetylacetonato 1,5-cyclooctadiene silver(I) (hfacAg<sup>I</sup>COD, **Cat. No. 348198**) is another early example of silver metal deposited by thermal ALD used by injection of the precursor dissolved in toluene, followed by a propanol wash. The publication does not describe the ALD process in detail, but suggests the propanol undergoes oxidative dehydrogenation (possibly catalysed by the silver metal) and that this is responsible for the reduction of the silver (**Figure 5**).



**Figure 5.** The oxidative dehydrogenation of propanol to produce a silver hydride moiety at a substrate surface.


The process was further investigated to determine the ALD process parameters and surface chemistry, highlighting several mechanistic points. First, it showed the cyclooctadiene is lost either in solution or during volatilization, emphasizing the surface moiety should be Ag<sup>I</sup>hfac (**Figure 5**). Second, the ALD temperature window was very small (between 121–130 °C) and any deviation adds a CVD element to the growth. Finally, the silver hydride surface was shown to improve the kinetics of monolayer formation when carried through to the introduction of the hfacAg<sup>I</sup>COD. This result was also predicted computationally to enhance coverage. Ultimately, the process results in a nominal growth rate of 0.16 Å/cycle when run at 125 °C, which is low considering that chemisorption might be partially aided by the existence of a surface hydride. This same precursor shows similar growth rates (0.20 Å/cycle) over a slightly larger temperature window (105–128 °C) when tertbutylhydrazine was used instead of propanol. The enhanced process characteristics are thought to be due to the thermal rearrangement of tertbutylhydrazine to hydrazine and isobutene, where hydrazine undergoes oxidative dehydrogenation much more readily than propanol:



A similar silver compound (Ag<sup>I</sup>(fod)PEt<sub>3</sub>) was shown to deposit silver metal using hydrogen plasma as the reducing agent. The process was viable over a small temperature range: at 120 °C, saturative growth was shown to be 0.3 Å/cycle, and at 140 °C, it was 0.4 Å/cycle. The same precursor was used with methylamine–borane (BH<sub>3</sub>(NMe<sub>2</sub>H)) as the reducing agent and

showed saturative growth between 104 °C and 130 °C with a rate of 0.3 Å/cycle. One benefit of using Ag<sup>I</sup>(fod)PEt<sub>3</sub> is the ease of its synthesis. Starting with AgO, the precursor can be made in 96% yield using standard Schlenk techniques.


Deposition of Ag<sup>I</sup>(fod)PEt<sub>3</sub> can also be achieved using a plasma generated from ammonia gas. Interestingly, when plasma ammonia is used, the authors report the growth rate increased significantly to 2.4 Å/cycle and the oxygen impurity decreased to 2%. Nitrogen impurity was shown to increase from 2% to 7% when using hydrogen plasma with the use of plasma ammonia. The authors attributed both the improved growth rate and nitrogen impurity to a longer lived amine defect on the silver surface that allows better nucleation of the silver precursor on the surface. In contrast, the hydrogen-terminated silver likely has a shorter lifetime, and so nucleation of Ag<sup>I</sup>(fod)PEt<sub>3</sub> is not as kinetically viable, resulting in a low growth rate (**Figure 6**).



**Figure 6.** The nucleation of Ag<sup>I</sup>(fod)PEt<sub>3</sub> was found to occur more readily with an amine-terminated surface than with a clean silver metal surface.

## Gold

Gold metal processes by ALD are very rare — only two have been published to date. The first used Me<sub>3</sub>Au<sup>III</sup>PMe<sub>3</sub> (**Figure 7**) and oxygen plasma, followed by water. This three step process has a growth rate of 0.5 Å/cycle and a narrow ALD window of 120–130 °C, with CVD occurring at higher temperatures. Interestingly, this process produces a film containing gold and phosphorous in the absence of the water pulse, and requires the water to convert the P<sub>2</sub>O<sub>5</sub> formed by oxygen plasma exposure to H<sub>3</sub>PO<sub>4</sub>, which becomes a volatile by-product. The second gold ALD process uses Me<sub>2</sub>Au<sup>III</sup>(S<sub>2</sub>CNET<sub>2</sub>) (**Figure 7**) and ozone to produce gold metal with a growth rate of 1.1 Å/cycle. In this process, the ALD window is broader and reaches a higher temperature (120–200 °C), and the growth rate is significantly higher than the first process. Another significant difference in these processes is in the precursors: Me<sub>3</sub>Au<sup>III</sup>PMe<sub>3</sub> is a liquid at room temperature that can be synthesized with an overall 93% purified yield, while Me<sub>2</sub>Au<sup>III</sup>(S<sub>2</sub>CNET<sub>2</sub>) is a room temperature solid with an overall purified yield of ~10%. However, gold metal is now accessible from both plasma and thermal ALD processes, with reasonable growth rates. Although there is still much room for improvement, this indicates that gold can now be considered an ALD-accessible metal.



**Figure 7.** The two currently available processes for the deposition of gold metal by ALD.

## Conclusions

The three coinage metals (copper, silver and gold) are all accessible through multiple ALD processes. Copper processes are the most well-developed, due to the commercial use of copper metal interconnects on the nanometer scale, coupled with the straightforward synthetic chemistry of copper-containing inorganic and organometallic compounds. There are a large number of copper precursors with various ligands, and processes include thermal methods at a variety of temperatures as well as plasma processes that can be used to deposit copper metal reliably. Presently, the preferred copper precursor appears to be  $Cu^{II}(\text{DMAP})$ , due to its ease of synthesis, versatility of deposition processes, and the low impurity of the resulting films.

Silver and gold metal deposition processes are less mature than processes using copper. Silver metal has been deposited by a variety of  $\beta$ -diketonate precursors, with fluorinated ligands showing the best results. Several of the silver processes have been very well characterized, and this should allow additional improvement of silver-containing ALD precursors.

For gold, the only two currently available processes cover thermal and plasma, but development has really just begun. The first true gold metal process was introduced in 2016, with the second following in 2017. Further understanding of both processes is required, and further precursor development will certainly help address the specific challenges for deploying gold across a wider temperature window with increased possibilities for surface chemistry. The lessons learned about copper ALD can now be applied to deposition of silver and gold, and the possibilities are just beginning to be explored.

## References

- (1) Gordon, P. G.; Kurek, A.; Barry, S. T. *ECS JSS* **2015**, *4*, N3188–N3197.
- (2) Miikkulainen, V.; Leskelä, M.; Ritala, M.; Puurunen, R. L. *J. Appl. Phys.* **2013**, *113*, 021301.
- (3) Young, V. L.; Cox, D. F.; Davis, M. E. *Chem. Mater.* **1993**, *5*, 1701–1709.
- (4) Elliott, S. D.; Dey, G.; Maimaiti, Y. J. *Chem. Phys.* **2017**, *146*, 052822.
- (5) Lee, B. H.; Hwang J. K.; Nam J. W.; Lee S. U.; Kim J. T.; Koo S. M.; Baunemann, A.; Fischer, R. A.; Sung, M. M. *Angew. Chem. Int. Ed. Engl.* **2009**, *48*, 4536–4539.
- (6) Knisley, T. J.; Ariyasena, T. C.; Sajavaara, T.; Saly, M. J. Winter, C. H. *Chem. Mater.* **2011**, *23*, 4417–4419.
- (7) Kalutarage, L. C.; Clendenning, S. B.; Winter, C. H. *Chem. Mater.* **2014**, *26*, 3731–3738.
- (8) Niskanen, A.; Hatanpää, T.; Arstila, K.; Leskelä, M.; Ritala, M. *Chem. Vap. Deposition* **2007**, *13*, 408–413.
- (9) Masango, S. S.; Peng, L.; Marks, L. D.; Van Duyne, R. P.; Stair, P. C. *J. Phys. Chem. C* **2014**, *118*, 17655–17661.
- (10) Chalker, P. R.; Romani, S.; Marshall, P. A.; Rosseinsky, M. J.; Rushworth, S.; Williams, P. A. *Nanotechnology* **2010**, *21*, 1–7.
- (11) Golrokhi, Z.; Chalker, S.; Sutcliffe, C. J.; Potter, R. J. *Appl. Surf. Sci.* **2016**, *364*, 789–797.
- (12) Golrokhi, Z.; Marshall, P. A.; Romani, S.; Rushworth, S.; Chalker, P. R.; Potter, R. J. *Appl. Surf. Sci.* **2017**, *399*, 123–131.
- (13) Kariniemi, M.; Niinistö, J.; Hatanpää, T.; Kemell, M.; Sajavaara, T.; Ritala, M.; Leskelä, M. *Chem. Mater.* **2011**, *23*, 2901–2907.
- (14) Mäkelä, M.; Hatanpää, T.; Mizohata, K.; Meinander, K.; Niinistö, J.; Räisänen, J.; Ritala, M.; Leskelä, M. *Chem. Mater.* **2017**, *29*, 2040–2045.
- (15) Minjauw, M. M.; Solano, E.; Sree S. P.; Asapu, R.; Van Daele, M.; Ramachandran, R. K.; Heremans, G.; Verbruggen, S. W.; Lenaerts, S.; Martens, J. A.; Detavernier, C.; Dendooven, J. *Chem. Mater.* **2017**, *29*, 7114–7121.
- (16) Griffiths, M. B. E.; Pallister, P. J.; Mandia, D. J.; Barry, S. T. *Chem. Mater.* **2016**, *28*, 44–46.
- (17) Mäkelä, M.; Hatanpää, T.; Mizohata, K.; Räisänen, J.; Ritala, M.; Leskelä, M. *Chem. Mater.* **2017**, *29*, 6130–6136.

## Group 11 Precursors

For a complete list of available materials, visit [SigmaAldrich.com/solutiondeposition](http://SigmaAldrich.com/solutiondeposition).

### Copper

Name	Structure	Purity (%) Or Grade	Form	Cat. No.
Copper(I) acetate		97%	powder and chunks	<b>403342-1G</b> <b>403342-10G</b>
Copper(II) acetate		99.99% trace metals basis	powder	<b>517453-5G</b> <b>517453-25G</b>
		98%	powder or crystals	<b>326755-25G</b> <b>326755-100G</b>
Copper(II) acetate monohydrate		99.99% trace metals basis	powder or crystals	<b>229601-10G</b> <b>229601-50G</b>
Copper(II) acetate hydrate		98%	crystals	<b>341746-100G</b> <b>341746-500G</b> <b>341746-2.5KG</b>
Copper(II) acetylacetonate		≥99.9% trace metals basis	powder	<b>514365-10G</b> <b>514365-50G</b>
		97%	powder	<b>C87851-25G</b> <b>C87851-100G</b> <b>C87851-500G</b>
Copper(II) 2-ethylhexanoate		-	powder or chunks	<b>337323-5G</b> <b>337323-25G</b>
Copper(II) hexafluoroacetylacetonate hydrate		-	powder or chunks	<b>335193-5G</b>
Copper(II) trifluoroacetate hydrate		-	powder or chunks	<b>311189-5G</b>
Copper(II) trifluoroacetylacetonate		97%	solid	<b>101826-5G</b>

### Silver

Name	Structure	Purity (%) Or Grade	Form	Cat. No.
(1,5-Cyclooctadiene) (hexafluoroacetylacetonato)silver(I)		99%	powder	<b>348198-5G</b>
Silver acetate		99.99% trace metals basis	powder or crystals	<b>204374-10G</b> <b>204374-50G</b>
Silver acetylacetonate		98%	solid	<b>323489-1G</b> <b>323489-5G</b>
Silver benzoate		99%	powder or crystals	<b>227277-10G</b> <b>227277-50G</b>
Silver trifluoroacetate		≥99.99% trace metals basis	liquid	<b>482307-1G</b> <b>482307-2.5G</b> <b>482307-10G</b>
		98%	powder	<b>T62405-5G</b> <b>T62405-25G</b>

### Gold

Name	Structure	Purity (%) Or Grade	Form	Cat. No.
Gold(I) cyanide	AuCN	99.9% trace metals basis	powder	<b>254088-1G</b>




# Few Monolayer Atomic Layer Deposition (ALD) on Surfaces and Interfaces for Energy Applications



Peifu Cheng, Parag Banerjee\*

Department of Mechanical Engineering and Materials Science  
One Brookings Drive, Washington University in St. Louis, St. Louis, MO 63130  
\*Email: parag.banerjee@wustl.edu



## Introduction

Atomic layer deposition (ALD) exploits a sequential set of two or more surface reactions repeated cyclically to achieve thin film growth.<sup>1</sup> These surface reactions occur between a chemisorbed species attached to a substrate and incoming precursor molecules from the gas-phase. Due to its precise deposition mechanism, ALD possesses many advantages over other thin film deposition processes. First, ALD offers extremely precise control at the angstrom level over the thickness of the deposited layer. Second, ultrahigh-aspect-ratio structures, including nanoporous solids<sup>2</sup> and three-dimensional (3D) hierarchical structures,<sup>3</sup> can be conformally and uniformly covered by ALD without thickness or composition gradients. Third, the chemical composition of films can be controlled by tailoring ALD cycles in a programmable manner.<sup>4,5</sup> Finally, ALD enables deposition at lower temperatures,<sup>6</sup> making it possible to use with polymer-based<sup>7</sup> and biological substrates that are incompatible with higher temperatures.<sup>8</sup>

It is not widely appreciated that only a few monolayers of material added to a surface or interface can significantly impact performance, quality, and reliability. Since ALD allows atomic scale resolution in deposition of monolayers, new scientific mechanisms and phenomena can be investigated by characterizing these surfaces. While there are many excellent review articles on applied ALD research,<sup>9-17</sup> none focus on the impact of few monolayer ALD films on materials and devices. Therefore, this review highlights research using few monolayer ALD to enhance material and device quality, performance, and reliability, specifically in the fields of photovoltaics, batteries, and catalysis.


We define “few monolayers” as less than 10 layers or a thickness of less than 2–3 nm. This demarcation of a ‘few’ is based on research reviews, which show impressive improvements in performance and reliability within this thickness range, and in most instances actual performance suppression by a thicker film. This definition also includes descriptions such as ultrathin, thin, compact, and few layer.

## Photovoltaics

ALD has been used to suppress charge recombination at interfaces in photovoltaic devices such as hybrid organic-inorganic perovskite solar cells, dye-sensitized solar cells (DSSCs), and quantum dot solar cells (QDSCs). Also, with the advent of hybrid organic-inorganic perovskite solar cells, ALD has been effectively used to improve the stability of these air sensitive compounds.

TiO<sub>2</sub> (up to four nm) was deposited by ALD on spin-coated mesoporous nanoparticulate TiO<sub>2</sub> films as blocking layers to make high efficiency solid-state perovskite-based (CH<sub>3</sub>NH<sub>3</sub>PbI<sub>3</sub>) solar cells.<sup>18</sup> As shown in **Figure 1A**, a two nm conformal TiO<sub>2</sub> ALD overlayer was sufficient to block the electron recombination, leading to a large open-circuit potential of 969.4 mV (reference potential 836.6 mV) and the high photon-to-current efficiency (PCE) of 11.5%.

Few layer ALD has been employed to enhance the PCE of DSSCs.<sup>19-22</sup> An ultrathin TiO<sub>2</sub> layer (10 cycles, 0.3–0.6 nm) was successfully coated on the surface of submicrometer-sized aggregates of ZnO nanocrystallites.<sup>19</sup> As shown in **Figure 1B**, both the open-circuit voltage and the fill factor were increased



**Figure 1.** A) Two nm ALD TiO<sub>2</sub> on mesoporous TiO<sub>2</sub> is sufficient to block electron hole recombination in solid-state perovskite-based solar cells. Reproduced with permission from reference 18, copyright 2014 Wiley. B) 0.6 nm of ALD TiO<sub>2</sub> on nanocrystalline ZnO suppresses charge recombination in DSSCs. Reproduced with permission from reference 19, copyright 2010 Wiley. C) 0.3 nm Al<sub>2</sub>O<sub>3</sub> improves ambient stability of hybrid organic-inorganic perovskite solar cells to up to 24 days. Reproduced with permission from reference 24, copyright 2015 Royal Society of Chemistry.

by the introduction of a TiO<sub>2</sub> layer, due to suppression of surface charge recombination, without impairing the photocurrent density. This resulted in a more than 20% enhancement of PCE from 5.2% to 6.3%.

Conformal TiO<sub>2</sub> interfacial layers with a thickness of 1.5, 2.2, and 3.2 nm were coated on CdS quantum dots (QDs) and N719 dyes to enhance the performance of quantum dot and dye co-sensitized solar cells.<sup>23</sup> The introduction of TiO<sub>2</sub> interfacial layers significantly enhanced the short-circuit current density and PCE by improving the stability of the QDs and reducing the electron recombination with the electrolyte. The highest PCE reported was 2.36%, a 41% improvement from the reference device, at the optimal TiO<sub>2</sub> thickness of 2.2 nm.

Ultrathin compact Al<sub>2</sub>O<sub>3</sub> layers were introduced by ALD (1-7 cycles) onto a layer of hole-transport material (HTM) to improve the stability of hybrid organic-inorganic perovskite solar cells in humid conditions.<sup>24</sup> The experimental results shown in **Figure 1C** indicate that a 3-cycle (~0.3 nm) ALD Al<sub>2</sub>O<sub>3</sub> layer greatly improved the ambient stability of organic-inorganic perovskite solar cells and retained 90% of the initial PCE after 24 days of storage in air. When the number of ALD cycles increased, the current decreased sharply, probably caused by a significant reduction in the quantum tunneling of the holes.

## Batteries


Batteries face limitations in energy storage density, charging rate, cycling stability, safety, and cost.<sup>17</sup> Few layer ALD can effectively address many of these issues. For example, nanosized LiCoO<sub>2</sub> powders were prepared by a molten salt method, and then coated directly using ALD (2 and 6 cycle) to make ultrathin conformal Al<sub>2</sub>O<sub>3</sub> films (**Figure 2A**).<sup>25</sup> The nano-LiCoO<sub>2</sub> electrodes coated with 2 ALD cycles of Al<sub>2</sub>O<sub>3</sub> (0.22 nm) delivered the highest discharge capacity of 133 mAh/g at a current rate of 1400 mA/g (7.8 C), corresponding to a 250% improvement in reversible capacity when compared to bare LiCoO<sub>2</sub> nanoparticles. The improved performance resulted from the suppression of cobalt dissolution from the LiCoO<sub>2</sub> by the

ALD Al<sub>2</sub>O<sub>3</sub> film. However, coating nano-LiCoO<sub>2</sub> electrodes with 6 ALD cycles of Al<sub>2</sub>O<sub>3</sub> decreased the rate performance, indicating that thick films resulted in poor Li<sup>+</sup> conductivity.

Nanostructured anodes coated by few layer ALD have also shown enhanced rate performance and cycling stability.<sup>26</sup> Porous Fe<sub>2</sub>O<sub>3</sub> nanorods anchored on nitrogen-doped graphene sheets (NGr) (**Figure 2B**) are one example. ALD was used to deposit two cycles of Al<sub>2</sub>O<sub>3</sub> film (<1 nm) on the Fe<sub>2</sub>O<sub>3</sub> nanorods anchored on NGr (NGr-M-I-2ALD). This anode displayed a first-cycle Coulombic efficiency of 89%, much higher than uncoated NGr-I-M (65%). Further, the stable capacity of the anodes (249 mAh/g) at 20 A/g was maintained for 200 cycles.

Few layer ALD was also used for cathode modification to improve the performance of a lithium-sulfur (Li-S) battery.<sup>27</sup> Al<sub>2</sub>O<sub>3</sub> was deposited by ALD on a graphene-sulfur (G-S) composite electrode. The specific capacity (646 mAh/g) of a G-S composite cathode with 10 ALD cycles of Al<sub>2</sub>O<sub>3</sub> film (~1 nm) after 100 charge/discharge cycles at 0.5 oC was 103% higher than that of bare G-S (**Figure 2C**). The coulombic efficiency, rate capability, and electrochemical stability of the G-S composite electrode were also significantly increased by ALD-Al<sub>2</sub>O<sub>3</sub> coating.

Other examples of the use of few monolayer ALD in lithium ion batteries include direct deposition of Al<sub>2</sub>O<sub>3</sub> on natural graphite particles (**Figure 2D**).<sup>28</sup> These electrodes exhibit improved capacity retention when compared to bare electrodes because the insulating Al<sub>2</sub>O<sub>3</sub> layer blocks electron conduction paths between natural graphite particles and the current collector. Ultrathin Al<sub>2</sub>O<sub>3</sub> layers were coated directly on Li metal foils without air exposure to produce anodes for Li symmetric cells. The Li metal anodes treated with 20 ALD cycles of Al<sub>2</sub>O<sub>3</sub> (~2 nm) were used for up to 1259 cycles before failure, double the lifetime of non-treated anodes (control sample = 711 cycles) (**Figure 2E**).<sup>29</sup> Few layer ALD also helped prevent dendrite formation and improve the homogeneity of current distribution along the Li electrode/solid electrolyte interphase (SEI) interface.



**Figure 2.** A)  $0.2 \text{ nm Al}_2\text{O}_3$  coating on nanosized  $\text{LiCoO}_2$  leads to increase capacity and stability. Reproduced with permission from reference 25, copyright 2011 American Chemical Society. B) High stability of  $0.2 \text{ nm Al}_2\text{O}_3$  coated on  $\text{Fe}_2\text{O}_3$  attached to N-doped graphene sheets (NGr-I-M-2ALD) maintained capacity for 200 cycles. Reproduced with permission from reference 26, copyright 2014 Elsevier. C)  $1 \text{ nm Al}_2\text{O}_3$  deposited on graphene sulphur composite electrodes shows 103% high specific capacity compared to uncoated electrodes. Reproduced with permission from reference 27, copyright 2014 Royal Society of Chemistry. D) Direct deposition of  $\text{Al}_2\text{O}_3$  on graphitic electrodes shows stable capacity behavior. Reproduced with permission from reference 28, copyright 2010 Wiley. E)  $2 \text{ nm Al}_2\text{O}_3$  on Li metal foils shows improved cyclability and prevents dendrite formation.

## Catalysis

The challenges in catalyst development are improving catalytic activity, selectivity, and stability with minimal byproducts.<sup>14</sup> Current synthetic methods, which include impregnation, ion exchange, and precipitation, are limited in facing these challenges because of the non-uniform distribution of surface sites and the relatively narrow set of starting materials available. Few monolayer ALD has been tested as a way to address these challenges.

Metal oxide catalysts with high catalytic activity and stability have been developed by few monolayer ALD.<sup>30</sup>  $\text{TiO}_2$  thin films ( $15 \text{ nm}$ ) were deposited by ALD on supported Ag nanoparticles, while  $\text{SiO}_2$  thin interlayers with thicknesses of  $2, 5, 10,$  and  $20 \text{ nm}$  were deposited by ALD to separate the  $\text{TiO}_2$  film from Ag nanoparticles. The photocatalytic activities of the samples and references were investigated by degradation of methylene blue (MB) in UV light illumination (Figure 3A).  $\text{TiO}_2/\text{SiO}_2/\text{Ag}$  architectures with a  $2 \text{ nm SiO}_2$  interlayer displayed the best photocatalytic performance due to blue-shifting in the localized


surface plasmon resonance (LSPR) and coupling energy from Ag to  $\text{TiO}_2$  bandgap. Thicker  $\text{SiO}_2$  interlayers decreased the photocatalytic activity because they reduced the penetration depth of the LSPR electromagnetic field into the  $\text{SiO}_2$  layer.

Few layer ALD was also employed for the photocatalytic conversion of  $\text{CO}_2$  to  $\text{CO}$ .<sup>31</sup> Dense  $\text{ZnO}$  islands with different thickness were coated via ALD onto large arrays of massively parallel, single crystalline, high-density  $\text{CuO}$  nanowires generated by thermal oxidation of high-purity copper foils. The  $\text{CO}$  yield first increased and then decreased with increasing ALD cycles, and  $\text{CuO}$  nanowires with eight ALD cycles of  $\text{ZnO}$  islands (1.4 nm) displayed the optimal performance (**Figure 3B**). This behavior occurred because 1) the favorable band alignment and epitaxial matching between ALD  $\text{ZnO}$  nanoislands and  $\text{CuO}$  nanowires led to recombination-free charge transfer, 2) the surface morphology allowed gas exposure to both  $\text{ZnO}$  and  $\text{CuO}$  surfaces, and 3) surface defects on exposed  $\text{ZnO}$  islands resulted in lower quantum yields due to the trapping of transferred electrons.

The catalytic activity, selectivity, and stability of supported metal catalysts were also improved by few layer ALD.<sup>32–34</sup> One of the most important catalysts is supported platinum. Platinum nanoparticles were deposited on  $\text{Al}_2\text{O}_3$  by ALD (one or five

cycles) with  $\text{H}_2$  or  $\text{O}_2$  as the second half-cycle precursor for oxidative dehydrogenation of propane. Using  $\text{H}_2$  as the second precursor resulted in smaller Pt particles with lower average coordination numbers than  $\text{O}_2$ . After one ALD cycle (1.1 nm, measured by TEM), Pt catalyst with  $\text{H}_2$  as the second precursor exhibited the highest  $\text{C}_3\text{H}_6$  selectivity (37% at 14% conversion at 400 °C), likely due to the small particle size (**Figure 3C**). Commercial Pt catalysts (3.9 nm, measured by TEM) had a selectivity of less than 1% at 9% conversion under the same reaction conditions.<sup>32</sup>

$\text{TiO}_2$  overcoated  $\text{Co}/\text{TiO}_2$  catalysts were also prepared to replace precious-metal catalysts, like Pt and Ru, for aqueous-phase hydrogenation (APH) reactions.<sup>35</sup> Generally, base metals (like cobalt) are inexpensive, abundant, and active for liquid phase reactions, but they have low stability caused by irreversible deactivation. **Figure 3D** shows that a thin ALD  $\text{TiO}_2$  coating prevented both leaching and sintering of Co particles during aqueous-phase reactions, due to  $\text{TiO}_2$  decoration of under-coordinated cobalt sites located at defects, corners, and edges. The regenerated  $\text{TiO}_2/\text{Co}/\text{TiO}_2$  catalysts (30 ALD cycles, 1.2 nm) displayed much higher catalytic reactivity and stability than fresh catalysts for APH of furfuryl alcohol because more Co active sites were exposed after regeneration.<sup>35</sup>



**Figure 3.** A) Fast decay in methylene blue from  $\text{TiO}_2/\text{SiO}_2/\text{Ag}$  with 2 nm  $\text{SiO}_2$  separation layer (green decay curve). Reproduced with permission from reference 30, copyright 2011 American Chemical Society. B)  $\text{ZnO}$  islands on  $\text{CuO}$  nanowire have optimal size (1.4 nm) for maximal  $\text{CO}_2 \rightarrow \text{CO}$  conversion. Reproduced with permission from reference 31, copyright 2015 American Chemical Society. C) Pt nanoparticle catalyst (1.1 nm) using  $\text{H}_2$  as reactant shows highest selectivity to oxidative dehydrogenation to propane at 400 °C. Reproduced with permission from reference 32, copyright 2015 American Chemical Society. D) 1.2 nm  $\text{TiO}_2$  on Co on  $\text{TiO}_2$  showed highest activity and stability for aqueous phase hydrogenation reactions. Reproduced with permission from reference 35, copyright 2014 Royal Society of Chemistry.



## Conclusions

In summary, ALD is a singularly powerful technique that can be applied to tailor surfaces and interfaces using few monolayer films, which have a significant impact on heterogeneous material and device performance. The structure of the few monolayers can be precisely engineered, depending on how the precursor molecule interacts with the chemisorbed ligands on the substrate. If the deposition is layer-by-layer, conformal films are produced that are ideal for charge blocking or for providing surface passivation and protection. ALD films using these methods have been applied in a number of photovoltaics and batteries. On the other hand, if the ALD chemistry is driven by island nucleation followed by growth, then surfaces and interfaces can be decorated with nanocrystals in order to store charge and present chemically active facets useful for inducing chemical reactions. Their diversity in tailoring surfaces and interfaces is unique to ALD processes.

As ALD processes continue to develop, new precursor molecules are needed to provide greater flexibility for tailoring surfaces and interfaces. The interaction of precursor ligands with substrate functional groups controls nucleation of ALD films at cycle 'one'. However, it is only now becoming clear that the interaction of nearest-neighbor adsorbed precursor molecules can affect film composition, structure, and thus properties.<sup>36-43</sup> This raises the interesting question as to whether the design of precursor molecules with appropriate ligand chemistry can be exploited to deterministically control film nucleation and growth on any surface.

Regardless of the application, ALD at its ultimate thickness limit has heralded a new era of surface and interface engineering. The control of this process appears to be remarkably simple, and yet has hidden complexities that will continue to push the boundaries of discovery of new materials and concept devices, and allow even finer control over their properties and applications.

## Acknowledgements

The authors acknowledge financial and material support from EMD Performance Materials. Helpful discussions with Dr. Ravi Kanjolia, Dr. Charles Dezelah and Dr. Jacob Woodruff are acknowledged. Partial funding was provided by U.S. Army RDECOM Acquisition Grant W911NF-15-1-0178, Subgrant RSC15032. Partial funding was also provided through the U.S.-India Partnership to Advance Clean Energy-Research (PACE-R) for the Solar Energy Research Institute for India and the United States (SERIUS), funded by the U.S. Department of Energy (Office of Science, Office of Basic Energy Sciences, and Energy Efficiency and Renewable Energy, Solar Energy Technology Program, under subcontract DE-AC36-08GO28308 to the national Renewable Energy Laboratory, Golden, Colorado).

## References

- (1) George, S. M. *Chem. Rev.* **2010**, *110*, 111-131.
- (2) Banerjee, P.; Perez, I.; Henn-Lecordier, L.; Lee, S. B.; Rubloff, G. W. *Nat. Nanotechnol.* **2009**, *4*, 292-296.
- (3) Bielineski, A. R.; Boban, B.; He, Y.; Kazyak, E.; Lee, D. H.; Wang, C. M.; Tuteja, A.; Dasgupta, N. P. *ACS Nano* **2017**, *11*, 478-489.
- (4) Banerjee, P.; Lee, W.-J.; Bae, K.-R.; Lee, S. B.; Rubloff, G. W. *J. Appl. Phys.* **2010**, *108*, 043504.
- (5) Gao, Z.; Myung, Y.; Huang, X.; Kanjolia, R.; Park, J.; Mishra, R.; Banerjee, P. *Adv. Mater. Interf.* **2017**, *3*, 1600496.
- (6) Groner, M. D.; Fabreguette, F. H.; Elam, J. W.; George, S. M. *Chem. Mater.* **2004**, *16*, 639-645.
- (7) Parsons, G. N.; Atanasov, S. E.; Dandley, E. C.; Devine, C. K.; Gong, B.; Jur, J. S.; Lee, K.; Oldham, C. J.; Peng, Q.; Spagnola, J. C.; Williams, P. S. *Coord. Chem. Rev.* **2013**, *257*, 3323-3331.
- (8) Knez, M.; Kadri, A.; Wege, C.; Gösele, U.; Jeske, H.; Nielsch, K. *Nano Lett.* **2006**, *6*, 1172-1177.
- (9) Knez, M.; Niesch, K.; Niinisto, L. *Adv. Mater.* **2007**, *19*, 3425-3438.
- (10) Kim, H.; Lee, H.-B.-R.; Maeng, W. J. *Thin Solid Films* **2009**, *517*, 2563-2580.
- (11) Parsons, G. N.; George, S. M.; Knez, M. *MRS Bull.* **2011**, *36*, 865-871.
- (12) Elam, J. W.; Dasgupta, N. P.; Prinz, F. B. *MRS Bull.* **2011**, *36*, 899-906.
- (13) Johnson, R. W.; Hultqvist, A.; Bent, S. F. *Mater. Today* **2014**, *17*, 236-246.
- (14) O'Neill, B. J.; Jackson, D. H. K.; Lee, J.; Canlas, C.; Stair, P. C.; Marshall, C. L.; Elam, J. W.; Kuech, T. F.; Dumesic, J. A.; Huber, G. W. *ACS Catal.* **2015**, *5*, 1804-1825.
- (15) Dasgupta, N. P.; Lee, H.-B.-R.; Bent, S. F.; Weiss, P. S. *Chem. Mater.* **2016**, *28*, 1943-1947.
- (16) Dasgupta, N. P.; Li, L.; Sun, X. *Adv. Mater. Interfaces* **2016**, *3*, 1600914.
- (17) Meng, X.; Wang, X.; Geng, D.; Ozgit-Akgun, C.; Schneider, N.; Elam, J. W. *Mater. Horizons* **2017**, *4*, 133-154.
- (18) Chandiran, A. K.; Yella, A.; Mayer, M. T.; Gao, P.; Nazeeruddin, M. K.; Grätzel, M. *Adv. Mater.* **2014**, *26*, 4309-4312.
- (19) Park, K.; Zhang, Q.; Garcia, B. B.; Zhou, X.; Jeong, Y.-H.; Cao, G. *Adv. Mater.* **2010**, *22*, 2329-2332.
- (20) Shanmugam, M.; Baroughi, M. F.; Galipeau, D. *Thin Solid Films* **2010**, *518*, 2678-2682.
- (21) Chandiran, A. K.; Tetreault, N.; Humphry-Baker, R.; Kessler, F.; Baranoff, E.; Yi, C.; Nazeeruddin, M. K.; Grätzel, M. *Nano Lett.* **2012**, *12*, 3941-3947.
- (22) Alibabaei, L.; Farnum, B. H.; Kalanyan, B.; Brennaman, M. K.; Losego, M. D.; Parsons, G. N.; Meyer, T. J. *Nano Lett.* **2014**, *14*, 3255-3261.
- (23) Xiu, L.; Kehan, Y.; Ganhua, L.; Junhong, C.; Chris, Y. J. *Phys. D* **2013**, *46*, 024004.
- (24) Dong, X.; Fang, X.; Lv, M.; Lin, B.; Zhang, S.; Ding, J.; Yuan, N. J. *Mater. Chem. A* **2015**, *3*, 5360-5367.
- (25) Scott, I. D.; Jung, Y. S.; Cavanagh, A. S.; Yan, Y.; Dillon, A. C.; George, S. M.; Lee, S.-H. *Nano Lett.* **2011**, *11*, 414-418.
- (26) Hu, T.; Xie, M.; Zhong, J.; Sun, H.-T.; Sun, X.; Scott, S.; George, S. M.; Liu, C.-S.; Lian, J. *Carbon* **2014**, *76*, 141-147.
- (27) Yu, M.; Yuan, W.; Li, C.; Hong, J.-D.; Shi, G. J. *Mater. Chem. A* **2014**, *2*, 7360-7366.
- (28) Jung, Y. S.; Cavanagh, A. S.; Riley, L. A.; Kang, S.-H.; Dillon, A. C.; Groner, M. D.; George, S. M.; Lee, S.-H. *Adv. Mater.* **2010**, *22*, 2172-2176.
- (29) Kazyak, E.; Wood, K. N.; Dasgupta, N. P. *Chem. Mater.* **2015**, *27*, 6457-6462.
- (30) Kumar, M. K.; Krishnamoorthy, S.; Tan, L. K.; Chiam, S. Y.; Tripathy, S.; Gao, H. *ACS Catal.* **2011**, *1*, 300-308.
- (31) Wang, W.-N.; Wu, F.; Myung, Y.; Niedzwiedzki, D. M.; Im, H. S.; Park, J.; Banerjee, P.; Biswas, P. *ACS Appl. Mater. Interfaces* **2015**, *7*, 5685-5692.
- (32) Gould, T. D.; Lubers, A. M.; Corpuz, A. R.; Weimer, A. W.; Falconer, J. L.; Medlin, J. W. *ACS Catal.* **2015**, *5*, 1344-1352.
- (33) Yan, H.; Cheng, H.; Yi, H.; Lin, Y.; Yao, T.; Wang, C.; Li, J.; Wei, S.; Lu, J. J. *Am. Chem. Soc.* **2015**, *137*, 10484-10487.
- (34) Wang, H.; Wang, C.; Yan, H.; Yi, H.; Lu, J. J. *Catal.* **2015**, *324*, 59-68.
- (35) Lee, J.; Jackson, D. H. K.; Li, T.; Winans, R. E.; Dumesic, J. A.; Kuech, T. F.; Huber, G. W. *Energy Environ. Sci.* **2014**, *7*, 1657-1660.
- (36) Shi, L.; Gao, Z.; Liu, Z.; Myung, Y.; Banerjee, P. *Chem. Mater.* **2017**, *29*, 5458-5462.
- (37) Shirazi, M.; Elliott, S. D. *J. Comput. Chem.* **2014**, *35*, 244-259.
- (38) Tanskanen, J. T.; Hagglund, C.; Bent, S. F. *Chem. Mater.* **2014**, *26*, 2795.
- (39) Shirazi, M.; Elliott, S. D. *Nanoscale* **2015**, *7*, 6311-6318.
- (40) Weckman, T.; Laasonen, K. *Phys. Chem. Chem. Phys.* **2015**, *17*, 17322-17334.
- (41) Maimaiti, Y.; Elliott, S. D. *Chem. Mater.* **2016**, *28*, 6282-6295.
- (42) Goehry, C.; Schneider, N. J. *Phys. Chem. C* **2017**, *121*, 5871-5881.
- (43) Guerra-Nuñez, C.; Döbeli, M.; Michler, J.; Utke, I. *Chem. Mater.* **2017**, *29*, 8690-8703.



## Titania Nanomaterials for Support

For a complete list of available materials, visit [SigmaAldrich.com/titaniasupport](http://SigmaAldrich.com/titaniasupport).

Name	Particle Size	Purity	Form	Cat. No.
Titanium	particle size <100 nm	98.5% trace metals basis	dispersion nanoparticles	513415-5G
Titania paste, active opaque	avg. part. size 20 nm (active) avg. part. size ≤450 nm (scatter)	-	paste (cream)	791555-5G 791555-20G
Titania paste, reflector	avg. part. size 150–250 nm (scatter)	-	paste (white)	791539-5G 791539-20G
Titania paste, transparent	avg. part. size 20 nm (active)	-	paste (yellow)	791547-10G 791547-20G
Titanium dioxide	particle size 22–25 nm (BET) 22 nm	>95% (anatase (XRD))	paste (nanocrystalline colloid)	798495-25G
	particle size 18–20 nm (BET) 18 nm	-	paste (nanocrystalline colloid)	798509-25G
	particle size 22 nm & >150 nm (BET)	-	paste (nanocrystalline colloid)	798517-25G
	particle size 18–20 nm (BET)	>95% (anatase(XRD))	paste (nanocrystalline colloid)	798525-25G
Titanium(IV) oxide	primary particle size 21 nm (TEM)	≥99.5% trace metals basis	nanopowder	718467-100G
Titanium(IV) oxide, anatase	particle size <25 nm	99.7% trace metals basis	nanopowder	637254-50G 637254-100G 637254-500G
Titanium(IV) oxide, rutile	particle size <100 nm	99.5% trace metals basis	nanopowder, diam. × L ~10 × ~40 nm	637262-25G 637262-100G 637262-500G
Titanium(IV) oxide, mixture of rutile and anatase	particle size <50 nm (XRD) particle size <100 nm (BET)	99.5% trace metals basis	nanopowder	634662-25G 634662-100G
	particle size <250 nm (DLS) particle size ~21 nm (primary particle size of starting nanopowder)	99.9% trace metals basis	nanoparticle paste	700355-25G
	particle size ~21 nm (primary particle size of starting nanopowder) particle size <150 nm (volume distribution, DLS)	99.5% trace metals basis	dispersion nanoparticles	700347-25G 700347-100G
	particle size <100 nm particle size ~30 nm (primary particle size of starting nanopowder)	99.9% trace metals basis	dispersion nanoparticles	700339-100G




## Titania Precursors for Solar Cells

For a complete list of available materials, visit [SigmaAldrich.com/solutiondeposition](http://SigmaAldrich.com/solutiondeposition).

Name	Structure	Purity (%) Or Grade	Form	Cat. No.
Titanium diisopropoxide bis(acetylacetonate)		-	solution	325252-100ML 325252-500ML
Titanium(IV) isopropoxide		99.999% trace metals basis	liquid	377996-5ML 377996-25ML 377996-100ML
		97%	liquid	205273-4X25ML 205273-100ML 205273-500ML 205273-1L 205273-2L 205273-2.5L

## Lithium Precursors for Energy

For a complete list of available materials, visit [SigmaAldrich.com/solutiondeposition](http://SigmaAldrich.com/solutiondeposition).

Name	Structure	Purity (%) Or Grade	Form	Cat. No.
Lithium acetate		99.95% trace metals basis	powder or crystals	517992-100G
Lithium acetylacetonate		97%	solid	413046-25G
Lithium tert-butoxide		97%	powder and chunks	400173-25G 400173-100G



Name	Structure	Purity (%) Or Grade	Form	Cat. No.
Lithium <i>tert</i> -butoxide solution		-	liquid	398209-50ML 398209-250ML
		-	liquid	398195-50ML 398195-250ML
Lithium ethoxide	CH <sub>3</sub> CH <sub>2</sub> OLi	95%	powder and chunks	400203-5G
Lithium ethoxide solution	CH <sub>3</sub> CH <sub>2</sub> OLi	-	liquid	400254-100ML
Lithium isopropoxide		95%	powder and chunks	348937-25G
Lithium methoxide	CH <sub>3</sub> OLi	98%	powder	344370-25G 344370-100G
Lithium methoxide solution	CH <sub>3</sub> OLi	-	liquid	408794-100ML 408794-800ML
Lithium phenoxide solution		-	liquid	400378-100ML
Lithium thiophenolate solution		-	liquid	400270-100ML 400270-800ML

## Zinc Precursors for Catalysis

For a complete list of available materials, visit [SigmaAldrich.com/solutiondeposition](http://SigmaAldrich.com/solutiondeposition).

Name	Structure	Purity (%) Or Grade	Form	Cat. No.
Zinc acetate dihydrate		99.999% trace metals basis	powder or chunks	379786-5G 379786-25G
Zinc acetylacetonate hydrate		99.995% trace metals basis	solid	480991-5G 480991-25G
		-	solid	132306-100G 132306-500G
Zinc hexafluoroacetylacetonate dihydrate		98%	solid	417858-1G
Zinc trifluoroacetate hydrate		-	powder	394017-5G

# subscribe today

Don't miss another  
topically focused technical review.

It's **free** to sign up for a print or digital  
subscription of *Material Matters*™.

- Advances in cutting-edge materials
- Technical reviews on emerging technology from leading scientists
- Peer-recommended materials with application notes
- Product and service recommendations




# MERCK

To view the library of past issues  
or to subscribe, visit  
[SigmaAldrich.com/materialmatters](http://SigmaAldrich.com/materialmatters).

# BIOMATERIALS FOR MECHANOREGENERATION

**Dr. David Mooney**  
Named 2018 Materials Research Society  
Mid-Career Researcher



Mechanical cues regulate many aspects of biology, and we are developing materials capable of providing defined mechanical input to cells and tissues. Relevant signals include both resistance to cell intrinsic mechanical forces and application of extrinsic forces. These biomaterials can enhance the effectiveness of stem cell therapies, directly support and regenerate damaged tissues, and have led to bio-inspired new medical adhesives with unprecedented properties.

## About David Mooney

David Mooney is the Pinkas Family Professor of Bioengineering in the Harvard School of Engineering and Applied Sciences, and a Core Faculty Member of the Wyss Institute. His laboratory designs biomaterials to make cell and protein therapies effective and practical approaches to treat disease. He is a member of the National Academy of Engineering, the National Academy of Medicine, and the National Academy of Inventors. He has won numerous awards, including the Clemson Award from the SFB, MERIT award from the NIH, Distinguished Scientist Award from the IADR, Phi Beta Kappa Prize for Excellence in Undergraduate Teaching, and the Everett Mendelsohn Excellence in Mentoring Award from Harvard College. His inventions have been licensed by numerous companies, leading to commercialized products, and he is active on industrial scientific advisory boards.

The life science business of Merck operates as  
MilliporeSigma in the U.S. and Canada.