



Mobius® Breez Microbioreactor

A 2 mL Automated Perfusion
Cell Culture Platform





- Run more experiments in parallel for higher efficiency
- Reliable run-to-run performance, obtain high-quality perfusate on demand
- Freedom from hands-on work with automated control and monitoring
- Operational flexibility allows you to save on development time and costs
- Easy to set up and use, allowing you to focus on your goals

Working Volume	2 mL
# of Controllers	4 (independent)
Flow Rate	0.1 – 10 VVD
Sample Volume	50 µL (min.)
Closed Loop Control	pH, DO, T, density-OD

Do more with less and accelerate your path to process intensification

- ✓ Perfusion process development
- ✓ Cell line development
- ✓ Media screening and optimization



 **Ask your sales rep about a trial!**