



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

SHEEP IgG Purified Immunoglobulin

Product Number **I 8265**

Product Description

Purified Sheep IgG is isolated from pooled normal sheep serum by fractionation.

Technical grade IgG may be used as an economical alternative to the reagent grade immunoglobulins, in applications where high purity is not required.

Reagents

The product is supplied as a liquid in 0.01 M phosphate buffered saline, pH 7.2, containing 15 mM sodium azide as a preservative.

Precautions and Disclaimer

Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling practices.

Storage/Stability

Store, in aliquots, at -20°C . Repeated freezing and thawing of product is **not** recommended.

JWM/KMR 03/02

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.